

CR AND DR ACCESSORIES



**X-RAY PLATE
DISPOSABLE
SMALL PLASTIC
COVERS (PACK
OF 25) 2**

X-ray Plate
Protector. Small
Plastic Covers,
sold per pack of
25 bags.

Measurements:
330 X 360 mm

**CODE C-
GEN1001**



**X-RAY PLATE
DISPOSABLE
LARGE PLASTIC
COVERS (PACK
OF 25)**

X-ray Plate
Protector. Large
Plastic Covers,
sold per pack of
25 bags.

Measurements:
500 X 500 mm

**CODE C-
GEN1002**



**DR PROTECTIVE
TUNNEL**

The ATX DR
Protective Tunnel
is a protector
that can be used
to perform
podiatry x-ray
studies
eliminating
potential damage
that can occur.

This unit is
constructed to
protect a
35x43cm
(14x17inch) DR
or CR panel.

It has 6 rubber
feet to ensure a
strong footing
and provide grip
to all surfaces &
terrain.

**CODE E-
ACC1023.**



**X-RAY ON LIGHT

**CODE E-
ACC1027**