

CR AND DR ACCESSORIES



**X-RAY PLATE  
DISPOSABLE  
SMALL PLASTIC  
COVERS (PACK  
OF 25) 2**

X-ray Plate  
Protector. Small  
Plastic Covers,  
sold per pack of  
25 bags.

**CODE C-  
GEN1001**



**X-RAY PLATE  
DISPOSABLE  
LARGE PLASTIC  
COVERS (PACK  
OF 25)**

X-ray Plate  
Protector. Large  
Plastic Covers,  
sold per pack of  
25 bags.

**CODE C-  
GEN1002**



**DR PROTECTIVE  
TUNNEL**

The ATX DR  
Protective Tunnel  
is a protector  
that can be used  
to perform  
podiatry x-ray  
studies  
eliminating  
potential damage  
that can occur.

This unit is  
constructed to  
protect a  
35x43cm  
(14x17inch) DR  
or CR panel.

It has 6 rubber  
feet to ensure a  
strong footing  
and provide grip  
to all surfaces &  
terrain.

**CODE E-  
ACC1023.**



**X-RAY ON LIGHT  
\***

**CODE E-  
ACC1027**