

DENTAL



BITE BLOCK
COVER (300)

Disposable
Bite Block
Covers. Box of
300

CODE C-
DEN1001



APIXIA
PHOSPHOR
PLATE BARRIER
ENVELOPES

Made with
environmentally
friendly
materials
which avoid
cross
contaminati
on and
provide
protection
for the
patient.
Single use
barrier
protection
envelopes
for Apixia
PSP Plates.

Single use
barrier
protection
envelopes
for Apixia
PSP Plates.

- SIZE 0- 500 pieces
- SIZE 1- 500 pieces
- SIZE 2- 500 pieces
- SIZE 3- 500 pieces
- SIZE 4- 300 pieces
- SIZE 5- 200 pieces

CODE C-
DEN1002-7




APIXIA
PHOSPHOR
PLATES

- Comfortable & reusable phosphor plates.
- Imaging with phosphor plates offers all the flexibility and patient comfort of film.
- Phosphor plates are reusable a countless number of times to maximize your investment.

Four plate
sizes for
various
animal
structures.

- Size 0- 35x22 mm, 4 pieces
- Size 1- 40x24 mm, 4 pieces
- Size 2- 41x31



OPG CENTERING
POST

75x6mm OPG
Centering
Post.. Sold
individually

CODE C-
DEN1013

- Size 3-54x27mm, 4 pieces
- Size 4-76x57mm, 2 pieces

CODE C-DEN1009-12



OPG BITE BLOCK

OPG Bite Block to fit. Sold individually

CODE C-DEN1014



OPG EAR OLIVE

OPG Ear Olive, Pack of 50.

CODE C-DEN1015



COTTON DENTAL ROLL
PACK OF 40

Disposable Cotton Dental Rolls. Pack of 40

CODE C-DEN1017



STICKERS,
ADHESIVE
BITING TABS

CODE C-DEN1020


**SMART
POSITIONER FOR
PHOSPHOR
PLATE**


Flow's patented SMART system makes using phosphor plates just as easy, familiar, and convenient as film.

Simply load your phosphor plate into the SMART Cassette

Kit includes:

- 1 Anterior Aiming Arm
- 2 Anterior Aiming Rings
- 2 Posterior Aiming Arms
- 2 Posterior Aiming Rings
- 8 Universal Bite Blocks
- 20 Size 1 Cassettes
- 50 Size 2 Cassettes

**CODE C-
DEN1024**


**RAPID
PARALLELING
KIT WITH BITE
WING**

RAPiD is a right angle positioning device that assures correct film/cone alignment and prevents "cone cutting". It can be used for anterior, posterior, and bite-wing x-ray exams using any brand and size (3, 2, 1, and 0) of dental film.

**CODE C-
DEN1025**


**ALL-BITE
SENSOR HOLDER**

**CODE C-
DEN1026**


**UNI-VERSE-ALL
SENSOR
POSITIONER KIT**

Uni-Verse-All's unique sensor holders snap in and snap out. Use to take Anterior or Posterior/Bite Wing images. Each sensor holder adjusts "on the fly" to hold every size sensor, horizontally or vertically and are fully autoclavable.

**CODE C-
DEN1027**