

DENTAL



BITE BLOCK COVER (300)

Disposable Bite Block Covers. Box of 300

CODE C-DEN1001



APIXIA PHOSPHOR PLATE BARRIER ENVELOPES

Made with environmentally friendly materials which avoid cross contamination and provide protection for the patient. Single use barrier protection envelopes for Apixia PSP Plates.

Single use barrier protection envelopes for Apixia PSP Plates.

- SIZE 0- 500 pieces
- SIZE 1- 500 pieces
- SIZE 2- 500 pieces
- SIZE 3- 500 pieces
- SIZE 4- 300 pieces
- SIZE 5- 200 pieces

CODE C-DEN1002-7




APIXIA PHOSPHOR PLATES

- Comfortable & reusable phosphor plates.
- Imaging with phosphor plates offers all the flexibility and patient comfort of film.
- Phosphor plates are reusable a countless number of times to maximize your investment.

Four plate sizes for various animal structures.

- Size 0- 35x22 mm, 4 pieces
- Size 1- 40x24 mm, 4 pieces
- Size 2- 41x31



OPG CENTERING POST

75x6mm OPG Centering Post.. Sold individually

CODE C-DEN1013

- Size 3-54x27mm, 4 pieces
- Size 4-76x57mm, 2 pieces

CODE C-DEN1009-12



OPG BITE BLOCK

OPG Bite Block to fit. Sold individually

CODE C-DEN1014



OPG EAR OLIVE

OPG Ear Olive, Pack of 50.

CODE C-DEN1015



COTTON DENTAL ROLL
PACK OF 40

Disposable Cotton Dental Rolls. Pack of 40

CODE C-DEN1017



STICKERS,
ADHESIVE
BITING TABS

CODE C-DEN1020



SMART POSITIONER FOR PHOSPHOR PLATE


Flow's patented SMART system makes using phosphor plates just as easy, familiar, and convenient as film.

Simply load your phosphor plate into the SMART Cassette

Kit includes:

- 1 Anterior Aiming Arm
- 2 Anterior Aiming Rings
- 2 Posterior Aiming Arms
- 2 Posterior Aiming Rings
- 8 Universal Bite Blocks
- 20 Size 1 Cassettes
- 50 Size 2 Cassettes

CODE C-
DEN1024



RAPID PARALLELING KIT WITH BITE WING

RAPiD is a right angle positioning device that assures correct film/cone alignment and prevents "cone cutting". It can be used for anterior, posterior, and bite-wing x-ray exams using any brand and size (3, 2, 1, and 0) of dental film.

CODE C-
DEN1025



ALL-BITE SENSOR HOLDER

CODE C-
DEN1026



UNI-VERSE-ALL SENSOR POSITIONER KIT

Uni-Verse-All's unique sensor holders snap in and snap out. Use to take Anterior or Posterior/Bite Wing images. Each sensor holder adjusts "on the fly" to hold every size sensor, horizontally or vertically and are fully autoclavable.

CODE C-
DEN1027