

DENTAL POSITIONING AIDS



**STICKEMS,
ADHESIVE
BITEWING TABS**

**CODE C-
DEN1020**



**SMART
POSITIONER FOR
PHOSPHOR
PLATE**


Flow's patented SMART system makes using phosphor plates just as easy, familiar, and convenient as film.

Simply load your phosphor plate into the SMART Cassette

Kit includes:

- 1 Anterior Aiming Arm
- 2 Anterior Aiming Rings
- 2 Posterior Aiming Arms
- 2 Posterior Aiming Rings
- 8 Universal Bite Blocks
- 20 Size 1 Cassettes
- 50 Size 2 Cassettes

**CODE C-
DEN1024**



**RAPID
PARALLELING
KIT WITH BITE
WING**

RAPiD is a right angle positioning device that assures correct film/cone alignment and prevents "cone cutting". It can be used for anterior, posterior, and bite-wing x-ray exams using any brand and size (3, 2, 1, and 0) of dental film.

**CODE C-
DEN1025**



**ALL-BITE
SENSOR HOLDER**

**CODE C-
DEN1026**

**UNI-VERSE-ALL
SENSOR
POSITIONER KIT**

Uni-Verse-All's unique sensor holders snap in and snap out. Use to take Anterior or Posterior/Bite Wing images. Each sensor holder adjusts "on the fly" to hold every size sensor, horizontally or vertically and are fully autoclavable.

CODE C-
DEN1027