

DENTAL POSITIONING AIDS



STICKEMS, ADHESIVE BITEWING TABS

CODE C-
DEN1020



SMART POSITIONER FOR PHOSPHOR PLATE


Flow's patented SMART system makes using phosphor plates just as easy, familiar, and convenient as film.

Simply load your phosphor plate into the SMART Cassette

Kit includes:

- 1 Anterior Aiming Arm
- 2 Anterior Aiming Rings
- 2 Posterior Aiming Arms
- 2 Posterior Aiming Rings
- 8 Universal Bite Blocks
- 20 Size 1 Cassettes
- 50 Size 2 Cassettes

CODE C-
DEN1024



RAPID PARALLELING KIT WITH BITE WING

RAPiD is a right angle positioning device that assures correct film/cone alignment and prevents "cone cutting". It can be used for anterior, posterior, and bite-wing x-ray exams using any brand and size (3, 2, 1, and 0) of dental film.

CODE C-
DEN1025



ALL-BITE SENSOR HOLDER

CODE C-
DEN1026

**UNI-VERSE-ALL
SENSOR
POSITIONER KIT**

Uni-Verse-All's unique sensor holders snap in and snap out. Use to take Anterior or Posterior/Bite Wing images. Each sensor holder adjusts "on the fly" to hold every size sensor, horizontally or vertically and are fully autoclavable.

CODE C-
DEN1027