

FOAMS


Imaging Positioning, Patient Positioning, Support Foams,
Manufactured by ATX



POSITIONING
FOAM CIRCLE
200X50MM

Circle Foam:
Positioning
Foam circle
with a 200mm
diameter and
50mm
deep. (Raw or
covered in
Dartex or
Vinyl).


CODE c-FOAM02



POSITIONING
FOAM SMALL
WEDGE
210X75X75MM

Small Foam
Wedge: With a
triangular face
at each end
75x75x10mm.
(Raw or
covered in
Dartex or
Vinyl).

CODE c-FOAM03



POSITIONING
FOAM 15°
WEDGE
250X250X50MM

Foam 15°
Wedge: with
a square
250x250x50m
m base.


CODE c-FOAM04



POSITIONING
FOAM 90°
WEDGE
250X180X180M
M

Foam 90°
Wedge:
With a
triangular face
at each end
180x180x260
mm.

CODE c-FOAM05



POSITIONING
FOAM LONG 90°
WEDGE
900X180X180M

Foam 90° wedge with a triangular face at each end 900x180x180 mm.

CODE c-FOAM06



POSITIONING
FOAM THIN
BLOCK
250X200X50MM

Rectangular
Foam Thin Block,
250x200mm,
50mm in depth.

CODE c-FOAM07



POSITIONING
FOAM 90°
WEDGE
600X180X180M

600mm Long
Foam 90° wedge with a 2 triangular face s at the ends 180x180x250 mm.

CODE c-FOAM08



POSITIONING
FOAM 60°
WEDGE
900X250X180M

Foam long 60° wedge 900x250 x180mm


CODE c-FOAM09



POSITIONING
FOAM
ARTHROGRAM
PIECE
450X295X300M

Positioning Foam
Arthrogram,
450x295x300mm.
Raw

CODE c-FOAM10



POSITIONING
SET OF MIXED
FOAMS (SET 1)

Foam
wedges and
blocks are a
key
component
to creating
high-
quality radi
ographs.

- 2x
FOAM0
2 - Circle:
100x50
mm
- 2x
FOAM0
3 - Small
Wedge
210x75
x75mm
- 2x
FOAM0
4 - 15°
Wedge
:




POSITIONING
SET OF MIXED
FOAMS (SET 2)

SET OF
MIXED
FOAMS

Variety of
shapes and
sizes
available, or
design your
own.

- 2x
FOAM0
2 - Circle
- 2x
FOAM0
3 - Small
Wedge
- 2x
FOAM0
4 - 15°
Wedge
- 1x
FOAM0



POSITIONING
FOAM, FINGER
STEPPED,
115X55X30MM

Stepped Finger
Sponge RAW,
115x55x30mm

CODE c-FOAM13

250x250x50mm
• 2x
FOAM0
5 - 90°
Wedge
:

250x180x180mm
• 2x
FOAM0
7 -
Thin
Block:


250x200x50mm
• 2x
FOAM1
8 -
Block:

200x100x100mm

CODE C-
FOAM11-1

5 - 90°
Wedge
• 1x
FOAM0
6 -
Long 9
0°
Wedge
• 2x
FOAM0
7 -
Thin
Block

CODE C-
FOAM11-2




POSITIONING
FOAM, T.V.
POSITIONING
WEDGE,
510X450X100M
M

TV positioning
Wedge -
510x450x100mm

FOAM Options:

- Raw
- or Covered in
- Dartex
 - Vinyl

CODE c-FOAM14



POSITIONING "V"
CUT FOAM
BLOCK
230X240X150M
M

V Cut Foam Block
RAW -
230x240x150mm
.

CODE c-FOAM15



POSITIONING
FOAM,
MATTRESS
CUSTOM AS PER
DRAWING

Custom Mattress
as per request

CODE c-
FOAM16R



POSITIONING
FOAM, CUSTOM
FOAM PIECE

Custom Foam
Piece: Ordered,
as per request

CODE c-FOAM17



**POSITIONING
FOAM BLOCK
200X100X100M**

Foam block with
200x100x100mm
Square faces
at each end.

Foam Options:

- Raw
- Covered in Dartex,
- Covered in Vinyl,
- Covered in Dartex Sealed.

CODE c-FOAM18




**POSITIONING
FOAM BLOCK
180X180X600M**

Foam positioning
block with a
180x180mm
square face at
each end.

Foam Options:

- Raw
- Dartex
- Vinyl
- Sealed

CODE c-FOAM19




**POSITIONING
FOAM SQUARE
BLOCK
250X250X50MM**

Foam Square
block

Foam Options:

- Raw
- Dartex
- Vinyl
- Sealed

CODE c-FOAM20



**POSITIONING
FOAM TORSO
BLOCK
400X300X120M**

Foam Torso block
400x300x90mm
deep.


FOAM Options:

- Raw

or Covered in

- Dartex,
- Vinyl,

CODE C-FOAM21



**POSITIONING
SMALL ANIMAL
CRADLE**

Radiolucent
positioning
cradle for small,
medium and
large animals.

- Small:
300mm long
x 250mm
x 120mm
high.
- Medium:
400mm long
x 300mm
x 120mm
high.



**POSITIONING
FOAM, INFANT
IMMOBILIZER
COVERED WITH
DARTEX**

Infant
Immobilizer
covered in
DARTEX

1 x foam base
(450x400x50
mm) with 2 x
adjustable
bolsters
(450x100x90
mm)

fastened and
adjusted by
velcro
and covered



**POSITIONING
FOAM, HEAD
IMMOBILISER
250X270X180M**

Head Immobiliser
Foam:

is a U shaped
foam designed
to support the
head of a
patient during
X-ray, CT and
MRI
studies. 250x2
70x180mm

CODE c-FOAM61



**POSITIONING
FOAM,
FOREHEAD
IMMOBILISER
120X270X250M**

Forehead
Immobiliser:

Is a U
shaped
designed to
support the
forehead of a
patient
during X-ray,
CT and MRI
studies. Can
be covered
in Dartex or
vinyl. 120x2
70x250mm.



120mm
high
• Large:
500mm
long x
350mm
x
150mm
high

CODE C-FOAM24

with Dartex

CODE c-
FOAM60C

CODE C-FOAM62



RADIOGRAPHY
CALIPER

Easy to use
caliper for
anatomy
thickness
measurements.

Set you cm
thickness
generator factors
for optimum
plate dose.

CODE E-
ACC1033