

## X-RAY FULL SOLUTIONS



### AV CHOICE TABLE - FOR PORTABLE X- RAY

#### X-ray Table with Integrated Tubestand

- Solid weld constructed table with urine traps
- Integrated travelling tubestand
- Electric locks
- Foot activated exposure switch 3 x - one at each end and down the length of the table
- Bucky grid cabinet and tray

[Please contact us for pricing](#)

**CODE E-  
GEN1036**



### LONGTAIL WALL STAND

The VS50 wall stand has an amazing range of motion for the active imaging area that extends from 1.5 inches (3.8 cm) to 81 inches (205.7 cm) above the floor.

It enables operators to execute routine upright exams on any portion of the Longtail DR as well as lengthy long-length studies with a single exposure, without the risk of motion or stitching difficulties. For unsurpassed confidence and quality, the VS50 delivers a solid and secure stand in a fixed spot.

**CODE E-  
GEN1039**



### RAYMASTER FULL SPINE INTEGRATED DIGITAL X-RAY CHIROPRACTIC SOLUTION

The **Raymaster Full Spine Integrated Digital X-Ray Chiropractic solution** is an economical and small footprint radiographic system designed specifically for chiropractic care. It combines flexibility, efficiency and affordability. The system will instantly improve patient workflow by eliminating the use of cassettes and developing conventional x-ray film. Instead, our Universal Chiropractic Solution are powered by the intuitive ATX Longtail full spine DR system giving technologists access to high resolution digital imaging and processing capabilities. The Longtail is the world's first monolithic full-spine flat-panel DR detector. No more long scan wait, no image stitching, no image gaps, and no plate overlap. The acquisition time is sub one second in most cases - so it will substantially improve efficiency.

#### Key benefits:



### OCTAVE INTEGRATED DIGITAL X-RAY SOLUTION

OCTAVE VR series is next-generation digital radiography systems for veterinary applications. Equipped with cutting-edge technology that achieves high-image quality, improved system usability, and compact design, VR10X-40K delivers a fast and convenient work experience at your busy hospital environment for an affordable total investment. It is also an intuitive and ergonomic system, designed specifically for veterinary applications.

#### Features:

- Small Footprint: Octave is compactly designed to allow installation in smaller locations.
- Touch Screen Console: Enjoy the fast and simple preparation with a simple touch. You can control the console



- Mounts floor-to-wall or floor-to-ceiling
- Tube rotation enables precision alignment of the tube for non-angled, angled and upright views.
- Convenient display gives indication of tube rotation
- CMP 200 X-ray generator series offers advanced high quality radiographic capabilities with unmatched performance and reliability resulting in the lowest cost of ownership over the product lifetime
- Longtail Digital Radiology. It is the first of its kind to produce long-bone and full-spine digital studies instantaneously with the highest geometric accuracy.

**Lead time:**  
Subject to stock availability, transit times and radiation licensing. Estimated shipping date to be provided once the purchase is confirmed

**CODE E-**  
GEN1103

- for adjusting the angle of collimator, setting information for SID, and generator parameter (kVp, mA, mAs).
- 4-Way Movable Table: Simple and smooth adjustment of the table location is achieved with electromagnetic brakes, allowing easier, more convenient, and safer positioning of the patient
- Foot switch for easy control: Foot switch offers complete hands-free control of X-ray exposure, table lock, and collimator light for your convenient usage.

**Lead time:**  
Subject to stock availability, transit times and radiation licensing. Estimated shipping date to be provided once the purchase is confirmed

**CODE E-**  
GEN1110



### **IRAY CANIS DENTAL CAMERA X RAY**

The IRay Portable Camera XRay is a portable intraoral X-ray device.

It is equipped with a 7-inch LCD touch screen and Android operating system. It can utilise an iRay sensor, and images can be transferred directly via WiFi or USB, allowing for simple point-and-shoot functionality. The machine is compact and lightweight, with a variety of network connection options.

**CODE E-  
GEN1117**



### **RAYMASTER INTEGRATED DIGITAL X-RAY CHIROPRACTIC SOLUTION**

The **Raymaster Integrated Digital X-Ray Chiropractic solution** is an economical and small footprint radiographic system designed specifically for chiropractic care. It combines flexibility, efficiency and affordability. The system will instantly improve patient workflow by eliminating the use of cassettes and developing conventional x-ray film. Instead, our Universal Chiropractic Solution are powered by the intuitive ATX DR system giving technologists access to high resolution digital imaging and processing capabilities. The acquisition time is sub one second in most cases - so it will substantially improve efficiency.

#### **Key benefits of the Universal Upright Raymaster include:**

- Mounts floor-to-wall or floor-to-ceiling
- Tube rotation enables precision alignment of the tube for non-angled, angled and upright views.
- Convenient display gives



- indication of tube rotation
- CMP 200 X-ray generator series offers advanced high quality radiographic capabilities with unmatched performance and reliability resulting in the lowest cost of ownership over the product lifetime.

**Lead time:**  
Subject to stock availability, transit times and radiation licensing. Estimated shipping date to be provided once the purchase is confirmed

**CODE** EQU-GEN-10058