

LEAD LETTERS



**BARIUM PASTE
20ML**

Barium Sulphate, radiopaque paste that provides excellent radiographic visualization of the true hoof wall. The tube contains 20mls of paste. The paste appears as an opaque (white line) on the image, when applied in a continuous bead from the coronary band to the ground surface. Due to the very thin depth of horn at the most proximal crown of the hoof wall the actual horn margin is rarely visible on even the softest penetrating views.

**CODE C-
GEN1010**



**ID MARKERS, L
UNMOUNTED
LEAD LETTER
9MM**

**Unmounted
Lead Letter
(9mm)**

CODE C-ID1002



**ID MARKERS, R
UNMOUNTED
LEAD LETTER
9MM**

**Unmounted
Lead Letter
(9mm)**

CODE C-ID1003



**ID MARKERS, L
MOUNTED LEAD
LETTER 10MM**

**L Mounted Lead
Letter (9mm)**

CODE C-ID1004



**ID MARKERS, R
MOUNTED LEAD
LETTER 9MM**

"R" Mounted
9mm lead letter .
CODE C-ID1005



**L METAL
STENCIL 20MM
CODE C-ID1006**



**R METAL
STENCIL 20MM
CODE C-
ID1006-1**



**SADDLE
MARKER LR AP
CODE C-ID1008**



**ID MARKERS, L
& R LEAD
LETTER WITH 2
INITIALS
MOUNTED
CODE C-ID1012**



**ID MARKER OR
STENCIL,
CUSTOM
MARKER
Custom Lead
Marker - made on
request
CODE C-ID1014**