

SCALING TOOLS



**BARIUM PASTE  
20ML**

Barium Sulphate, radiopaque paste that provides excellent radiographic visualization of the true hoof wall. The tube contains 20mls of paste. The paste appears as an opaque (white line) on the image, when applied in a continuous bead from the coronary band to the ground surface. Due to the very thin depth of horn at the most proximal crown of the hoof wall the actual horn margin is rarely visible on even the softest penetrating views.

**CODE C-  
GEN1010**



**AKUCAL IMAGE  
SCALING DEVICE  
( 25MM )  
(COMPLETE)**

25mm image  
scaling device

COMPLETE with  
all parts  
including-

- Suction cup base for Horizontal or Vertical use
- Flexible arm to position marker at anatomy level
- Aluminium tip avoids breakage from wear and tear
- Numerical size indicators visible in every image ( 25mm, 25.4mm, 30mm )
- Surgical Grade stainless steel marker
- Lever to engage or disengage suction

**CODE E-  
ACC1029-25**



**AKUCAL  
SCALING TOOL ( 25MM  
REPLACEMENT  
MARKER )**

25mm

Replacement  
stainless steel  
marker for Image  
Scaling Tool

EACH

**CODE E-  
ACC1030-25**




**AKUCAL IMAGE  
SCALING DEVICE  
( 30MM )  
(COMPLETE)**

25mm image  
scaling device

COMPLETE with  
all parts  
including-

- Suction cup base for Horizontal or Vertical use
- Flexible arm to position marker at anatomy level
- Aluminium tip avoids breakage from wear and tear
- Numerical size indicators visible in every image ( 25mm, 25.4mm, 30mm )
- Surgical Grade stainless steel marker
- Lever to engage or disengage suction

**CODE E-  
ACC1031-30**



**AKUCAL  
SCALING TOOL ( 30MM  
REPLACEMENT  
MARKER )**

30mm

Replacement  
stainless steel  
marker for Image  
Scaling Tool

EACH

**CODE E-**  
ACC1032-30