

SCALING TOOLS



**BARIUM PASTE
20ML**

Barium Sulphate, radiopaque paste that provides excellent radiographic visualization of the true hoof wall. The tube contains 20mls of paste. The paste appears as an opaque (white line) on the image, when applied in a continuous bead from the coronary band to the ground surface. Due to the very thin depth of horn at the most proximal crown of the hoof wall the actual horn margin is rarely visible on even the softest penetrating views.

**CODE C-
GEN1010**



**AKUCAL IMAGE
SCALING DEVICE
(25MM)
(COMPLETE)**

25mm image
scaling device

COMPLETE with
all parts
including-

- Suction cup base for Horizontal or Vertical use
- Flexible arm to position marker at anatomy level
- Aluminium tip avoids breakage from wear and tear
- Numerical size indicators visible in every image (25mm, 25.4mm, 30mm)
- Surgical Grade stainless steel marker
- Lever to engage or disengage suction

**CODE E-
ACC1029-25**



**AKUCAL
SCALING TOOL (25MM
REPLACEMENT
MARKER)**

25mm

Replacement
stainless steel
marker for Image
Scaling Tool

EACH

**CODE E-
ACC1030-25**




**AKUCAL IMAGE
SCALING DEVICE
(30MM)
(COMPLETE)**

25mm image
scaling device

COMPLETE with
all parts
including-

- Suction cup base for Horizontal or Vertical use
- Flexible arm to position marker at anatomy level
- Aluminium tip avoids breakage from wear and tear
- Numerical size indicators visible in every image (25mm, 25.4mm, 30mm)
- Surgical Grade stainless steel marker
- Lever to engage or disengage suction

**CODE E-
ACC1031-30**



**AKUCAL
SCALING TOOL (30MM
REPLACEMENT
MARKER)**

30mm

Replacement
stainless steel
marker for Image
Scaling Tool

EACH

CODE E-
ACC1032-30