

## SKIN MARKERS



### BEEKLEY X-SPOT® 101

#### BEEKLEY X-SPOT 101 BLUE

Ref: 101 Blue,  
(QTY: 150/box)

*X-SPOT 101: Identify and differentiate specific points of reference in conventional simulation treatment planning, chest, and trauma imaging.*

- 1.5mm lead-free pellet
- easy Pick-Up tab

**CODE C-**  
BEEK101



### BEEKLEY Y-SPOT® 102

#### BEEKLEY Y-SPOT 102 Green

Ref: 102 Green  
(QTY: 132/box)

*Identify and differentiate specific points of reference in conventional simulation treatment planning, chest, and trauma imaging.*

- 2.0mm lead-free pellet
- easy Pick-Up tab

**CODE C-**  
BEEK102



### BEEKLEY D-SPOT® ARROW 104

#### BEEKLEY D-SPOT® ARROW 104

Ref: 104 (QTY: 99/box)

- lead-free arrow with easy Pick-Up tab
- for marking nodes and outer canthus

**CODE C-**  
BEEK104



### BEEKLEY T-SPOT® LINE 111

#### BEEKLEY T-SPOT® LINE 111

Ref: 111, (QTY: 275cm cut to measure roll/box)

- .8mm diameter metallic line.
- for marking scars, nodes, field borders, and match lines

**CODE C-**  
BEEK111



**BEEKLEY L-SPOT® 112**

"L" Spot  
Radiopaque "L"  
Orientation  
marker  
65 per box  
**CODE C-  
BEEK112**



**BEEKLEY R-SPOT® 113**

"R" Spot  
Metallic "R"  
Orientation  
marker  
64 per box  
  
**CODE C-  
BEEK113**



**BEEKLEY  
GUIDELINES® CT  
BIOPSY GRID  
117**

*Guidelines® CT  
Biopsy  
Grid 117*  
**Ref: 117,  
(QTY:  
20/box)**

- 4" x 5" biops grid
- 7 non-metallic radiopaque lines
- open pane design for easy mark placement
- ideal for CT-guided biopsies of the neck and other small areas

**CODE C-  
BEEK117**



**BEEKLEY CT-SPOT® 119**

**BEEKLEY CT-SPOT® 119**  
**Ref: 119, (QTY:  
85/box)**

*Ideal for marking 3-point set-ups and isocenters, these unique markers image brightly during simulation without streaking or artifact.*

- 2.3mm non-metallic pellet
- for 3-points set-ups and isocenters

**CODE C-  
BEEK119**



**BEEKLEY CT-SPOT® 120**

Ideal for marking 3-point set-ups and isocenters, these unique markers image brightly during simulation without streaking or artifact.

**BEEKLEY CT-SPOT® 120, Ref: 120 (QTY: 85/box)**

- 4.0mm non-metallic pellet
- for 3-point set-ups and isocenters

**CODE C-BEEK120**



**BEEKLEY PINPOINT®128**

**PinPoint® for Image Registration 128**

**Ref: 128, (QTY: 30 / box)**

- ideal for image registration performed under CT and MRI
- conical design with a 1.27mm diameter center hole

**CODE C-BEEK128**



**MR-SPOT PACKETS™ 185**

**MR-SPOT Packets™ 185**

**Ref: 185, (QTY: 40/box)**

- 1.75cm Radiance® filled packet
- ideal for extremities, soft tissue masses, and breast MRI
- designed for optimum patient comfort

**CODE C-BEEK185**



**BEEKLEY LOCALIZERS™ 210**

**Localizer Localization Markers 210**

**Ref: 210, (QTY: 40/box)**

- 2.3mm low-density pellet
- ideal for identifying soft tissue masses on the neck

**CODE C-BEEK210**



**BEEKLEY LOCALIZERS™ 211**

**LOCALIZERS™ LOCALIZATION MARKERS 211**

**Ref: 211, (QTY: 40/box)**



**BEEKLEY POINTGUARDS® 215**

**PointGuards® 215**

**Ref: 215, (QTY: 75/box)**

- 1.5cm x 5.0cm linear adhesive covers for



**BEEKLEY GUIDELINES® CT BIOPSY GRID 217**

**Guidelines® CT Biopsy Grid 217**

**Ref: 217, (QTY: 15 / box)**

- 7" x



**BEEKLEY ACCUGRID®300 4.65"**

**AccuGrid® Specimen Localizing Grid System 300**

**Ref: 300, (QTY: 12 / box)**

- 4.65" radiolucent grid inside a



- 4.0mm low-density pellet
- Ideal for identifying areas of interest on abdominal cases, or with larger patients

CODE C-BEEK211

longer, more linear treatment planning marks

CODE C-BEEK215

- 7.5" biopsy grid
- 12 non-metallic radiopaque lines
- open pane design for easy mark placement
- ideal for CT-guided biopsies of abdomen and other large areas

CODE C-BEEK217

- sealable surgical bag
- for the safe imaging, transport, and processing of surgical breast specimens

CODE C-BEEK300



BEEKLEY ACCUGRID®302  
5.65"

AccuGrid® Specimen Localizing Grid System 302

Ref: 302, (QTY: 10 / box)

- 5.65" radiolucent grid inside a sealable surgical bag
- for the safe imaging, transport, and processing of surgical breast specimens
- expanded grid size for use with larger specimens



BEEKLEY CRADLES®309

Cradles® Needle Localization Wire Protectors 309

Ref: 309, (QTY: 12 / box)

- 12" x 3" translucent cover with foam adhesive border
- disposable, clear covers for protecting and stabilizing the needle localization wire

CODE C-BEEK309



BEEKLEY BELLA BLANKETS®310

Bella Blankets® protective coverlets 310

Provide sanitary protection for patients and help immobilize hard to position breasts with Bella Blankets. Improve Clinical Outcomes, Increase Patient Satisfaction, Save Money and Increase Revenues.

Ref: 310, (QTY: 50/box)

- 22.9cm x




BEEKLEY POINTGUARDS® 315

PointGuards® 315

Ref: 315, (QTY: 100 / box)

- 2.5cm circular adhesive covers
- for smaller treatment planning marks
- see-through cover for quicker and more precise placement

CODE C-BEEK315

|                                                                                                                                                                                                                                                                                                                                                                                       |                                                                                                                                                                                                                                                                                                                                                                                                                                           |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p>S</p> <p><b>CODE C-<br/>BEEK302</b></p>                                                                                                                                                                                                                                                                                                                                            |                                                                                                                                                                                                                                                                                                                                                                                                                                           | <p>27.9cm</p> <ul style="list-style-type: none"><li>intended for use on 18cm x 24cm recept or plate</li></ul> <p><b>CODE C-<br/>BEEK310</b></p>                                                                                                                                                                                                                                                                                                                                                                                |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|  <p><b>BEEKLEY<br/>POINTGUARDS®<br/>318</b></p> <p><b>PointGuards®<br/>318</b></p> <p><b>Ref: 318, (QTY:<br/>85/box)</b></p> <ul style="list-style-type: none"><li>2.5cm circular adhesive covers</li><li>crosshair design</li><li>for hard to find marks</li></ul> <p><b>CODE C-<br/>BEEK318</b></p> |  <p><b>BEEKLEY<br/>COLDSPOT® 500</b></p> <p><b>ColdSPOT®<br/>500</b></p> <p><b>Sold Individually</b></p> <ul style="list-style-type: none"><li>2.75" x 3.0" gel pack</li><li>self-adhesive side tabs</li><li>for targeted and effective cold treatment following minimally invasive breast biopsy procedures</li></ul> <p><b>CODE C-<br/>BEEK500</b></p> |  <p><b>BEEKLEY X-<br/>SPOT®601<br/>FLORALS®</b></p> <p><b>X-<br/>SPOT® 6<br/>01<br/>Florals®<br/>/ 101<br/>Blue</b></p> <p><b>Ref: 601<br/>Florals®<br/>Nipple<br/>Marker,<br/>(QTY:<br/>150/box)<br/>Available<br/>in blue<br/>(Ref:101)</b></p> <ul style="list-style-type: none"><li>1.5mm lead-free pellet</li><li>easy Pik-UP tab</li><li>compatible with both digital and analog equipment</li></ul> <p><b>CODE C-<br/>BEEK601</b></p> |  <p><b>BEEKLEY S-<br/>SPOT® 607<br/>FLORALS®</b></p> <p><b>S-SPOT® CT<br/>Simulation<br/>Marker 607<br/>Florals®</b></p> <p>Help define architectural distortion at the biopsy site by delineating the exact location of incision.</p> <p><b>Ref: 607, (QTY:<br/>363cm cut to<br/>measure<br/>roll/box)</b></p> <ul style="list-style-type: none"><li>33mm diameter radiopaque line</li><li>burnout-resistant brighter image</li><li>for marking scars and field borders</li></ul> <p><b>CODE C-<br/>BEEK607</b></p> |



BEEKLEY A-SPOT® 651 FLORALS®

A-SPOT® Light Image® palpable mass marker 651 Florals®

Ref: 651 Florals®, (QTY: 120 / box)

- uniform low-density see-through triangle
- highlights without obscuring
- compatible with both digital and analog equipment

CODE C-BEEK651



BEEKLEY O-SPOT® LIGHT IMAGE® 652 FLORALS

O-SPOT® Light Image® 652 Florals

Ref: 652 (QTY: 58 / box)

- 1mm uniform low-density see-through ring
- highlights without obscuring
- compatible with both digital and analog equipment

CODE C-BEEK652



BEEKLEY S-SPOT® 677 FLORALS®

S-SPOT® Light Image® Soft 'n Stretchy® 677 Florals®

Ref: 677 Florals®, (QTY: 280cm cut to measure roll / box)

- 1mm diameter uniform low-density see-through line
- pinch-free material
- compatible with both digital and analog equipment

CODE C-BEEK677



BEEKLEY V-SPOT® 703 FLORALS®

V-SPOT® Designed for Digital®

Ref: 703 Florals®, (QTY: 130 / box)

- 2.3mm non-metallic pellet
- easy Pik-UP tab
- designed for use with the most sensitive digital equipment

CODE C-BEEK703



**BEEKLEY  
TOMOSPOT®  
781**

**TOMOSPOT®  
FOR  
MARKING  
NIPPLES 781**

Distinctly identify and permanently document nipples, moles, scars, palpable masses, and areas of pain or concern on mammography images with BEEKLEYSPOT® skin markers for mammography. The result is clear communication between patient, technologist and radiologist regarding areas of interest and more precise interpretation of the mammography.

**Ref: 781,  
(QTY:  
92/box)**

- 2.3mm non-metallic pellet
- pinch-free material with adhesive-free center
- designed for use with 3D breast tomosynthesis equipment

**CODE C-  
BEEK781**



**BEEKLEY NEW  
TOMOSPOT®  
782**

**TOMOSPOT®  
FOR  
MARKING  
MOLES 782**

Provide certainty that densities on mammographic images highlighted by the O-shaped marker are moles and not areas of concern.

**Ref: 782,  
(QTY:  
54/box)**

- 1mm low-density see-through ring
- pinch-free material
- designed for use with 3D breast tomosynthesis equipment

**CODE C-  
BEEK782**



**BEEKLEY NEW  
TOMOSPOT®  
783**

**TOMOSPOT®  
FOR  
MARKING  
SCARS 783**

Distinctly identify and permanently document nipples, moles, scars, palpable masses, and areas of pain or concern on mammography images with BEEKLEYSPOT® skin markers for mammography. The result is clear communication between patient, technologist and radiologist regarding areas of interest and more precise interpretation of the mammography.

**Ref: 783,  
(QTY:  
273cm cut  
measure  
roll/box)**

- 1mm diameter low-density see-through line
- pinch-free material
- designed for use with 3D breast tomosynthesis equipment



**BEEKLEY NEW  
TOMOSPOT®  
784**

**TOMOSPOT®  
FOR  
MARKING  
PALPABLE  
MASSES 784**

Distinctly identify and permanently document nipples, moles, scars, palpable masses, and areas of pain or concern on mammography images with BEEKLEYSPOT® skin markers for mammography. The result is clear communication between patient, technologist and radiologist regarding areas of interest and more precise interpretation of the mammography.

**Ref: 784,  
(QTY:  
116/box)**

- low-density see-through triangle
- pinch-free material
- designed for use with 3D breast tomosynthesis equipment

CODE C-  
BEEK783

CODE C-  
BEEK784



**BEEKLEY NEW  
TOMOSPOT®  
785**

*TomoSPOT®  
for marking  
non-palpable  
areas of  
concern / pain  
marker 785*

*Distinctly  
identify and  
permanently  
document ni  
pples,  
moles, scars,  
palpable  
masses, and  
areas of pain  
or concern  
on  
mammograp  
hy images  
with  
BeekeleySPOT  
® skin  
markers for  
mammograp  
hy. The  
result is  
clear  
communicati  
on between  
patient,  
technologist  
and  
radiologist  
regarding  
areas of  
interest and  
more precise  
interpretatio  
n of the  
mammograp  
hy.*

**Ref: 785,  
(QTY:  
58/box)**

- low-density see-through square
- pinch-free material
- designed for use with 3D



**BEEKLEY N-  
SPOT® 790  
FLORALS®**

**N-SPOT®  
Designed  
for  
Digital®  
Soft n  
Stretchy  
®**

**Ref: 790  
Florals®,  
(QTY:  
94/box)**

- 2.3mm non-metallic pellet
- pinch-free material with adhesive-free center
- designed for use with the most sensitive digital equipment

**CODE C-  
BEEK790**



**BEEKLEY O-  
SPOT®791  
FLORALS®**

**O-SPOT®  
Designed  
for  
Digital®  
Soft n  
Stretchy  
® 791  
Florals®**

*Provide  
certainty  
that  
densities on  
mammograp  
hic images  
highlighted  
by the O-  
shaped  
marker are  
moles and  
not areas of  
concern.*

**Ref: 791  
(QTY:  
55/box)**

- 1mm uniform ultra-shear see-through ring
- pinch-free material
- designed for use with the most sensitive digital equipment

**CODE C-  
BEEK791**



**BARIUM PASTE  
20ML**

*Barium Sulphate,  
radiopaque paste  
that provides  
excellent  
radiographic  
visualization of  
the true hoof  
wall. The tube  
contains 20mls of  
paste. The paste  
appears as an  
opaque (white  
line) on the  
image, when  
applied in a  
continuous bead  
from the  
coronary band to  
the ground  
surface. Due to  
the very thin  
depth of horn at  
the most  
proximal crown  
of the hoof wall  
the actual horn  
margin is rarely  
visible on even  
the softest  
penetrating  
views.*

**CODE C-  
GEN1010**

*breast  
tomosy  
nthesis  
equip  
ment*

**CODE C-**  
BEEK785



**E-Z GRID FOR CT  
BIOPSY**

**E-Z Grid for CT  
biopsy. FAST  
FIND GRID.**

**CODE C-**  
MAR1012