


POSITIONING AIDS



**POSITIONING  
FOAM CIRCLE  
200X50MM**

**Circle Foam:**  
Positioning  
Foam circle  
with a 200mm  
diameter and  
50mm  
deep. (Raw or  
covered in  
Dartex or  
Vinyl).


**CODE** c-FOAM02



**POSITIONING  
FOAM SMALL  
WEDGE  
210X75X75MM**

**Small Foam  
Wedge:** With a  
triangular face  
at each end  
75x75x10mm.  
(Raw or  
covered in  
Dartex or  
Vinyl).

**CODE** c-FOAM03



**POSITIONING  
FOAM 15°  
WEDGE  
250X250X50MM**

**Foam 15°  
Wedge:** with  
a square  
250x250x50m  
m base.


**CODE** c-FOAM04



**POSITIONING  
FOAM 90°  
WEDGE  
250X180X180M  
M**

**Foam 90°  
Wedge:**  
With a  
triangular face  
at each end  
180x180x260  
mm.


**CODE** c-FOAM05



**POSITIONING  
FOAM LONG 90°  
WEDGE  
900X180X180M  
M**

**Foam 90° wedge**  
with a  
triangular face  
at each end  
900x180x180  
mm.

**CODE** c-FOAM06



**POSITIONING  
FOAM THIN  
BLOCK  
250X200X50MM**

**Rectangular  
Foam Thin Block,**  
250x200mm,  
50mm in depth.

**CODE** c-FOAM07



**POSITIONING  
FOAM 90°  
WEDGE  
600X180X180M  
M**

**600mm Long  
Foam 90° wedge**  
with a 2  
triangular face  
s at the ends  
180x180x250  
mm.

**CODE** c-FOAM08



**POSITIONING  
FOAM 60°  
WEDGE  
900X250X180M  
M**

**Foam long 60°  
wedge** 900x250  
x180mm


**CODE** c-FOAM09



POSITIONING  
FOAM  
ARTHROGRAM  
PIECE  
450X295X300M  
M

Positioning Foam  
Arthrogram,  
450x295x300mm.  
Raw

CODE c-FOAM10



POSITIONING  
SET OF MIXED  
FOAMS (SET 1)

Foam  
wedges and  
blocks are a  
key  
component  
to creating  
high-  
quality radi-  
ographs.

- 2x  
FOAM0  
2 - Circle:

100x50  
mm

- 2x  
FOAM0  
3 - Small  
Wedge

210x75  
x75mm

- 2x  
FOAM0  
4 - 15°  
Wedge

250x25  
0x50m  
m

- 2x  
FOAM0  
5 - 90°  
Wedge

250x18  
0x180m  
m

- 2x  
FOAM0  
7 - Thin  
Block:

250x20  
0x50m  
m

- 2x  
FOAM1  
8 - Block:

200x10  
0x100m  
m

CODE C-  
FOAM11-1



POSITIONING  
SET OF MIXED  
FOAMS (SET 2)

SET OF  
MIXED  
FOAMS

Variety of  
shapes and  
sizes  
available, or  
design your  
own.

- 2x  
FOAM0  
2 - Circle

- 2x  
FOAM0  
3 - Small  
Wedge


- 2x  
FOAM0  
4 - 15°  
Wedge

- 1x  
FOAM0  
5 - 90°  
Wedge

- 1x  
FOAM0  
6 - Long 9  
0°  
Wedge

- 2x  
FOAM0  
7 - Thin  
Block


CODE C-  
FOAM11-2



POSITIONING  
FOAM, FINGER  
STEPPED,  
115X55X30MM

Stepped Finger  
Sponge RAW,  
115x55x30mm

CODE c-FOAM13




**POSITIONING  
FOAM, T.V.  
POSITIONING  
WEDGE,  
510X450X100M  
M**

TV positioning  
Wedge -  
510x450x100mm

**FOAM Options:**

- Raw
- or Covered in
- Dartex
  - Vinyl
- CODE** c-FOAM14



**POSITIONING "V"  
CUT FOAM  
BLOCK  
230X240X150M  
M**

V Cut Foam Block  
RAW -  
230x240x150mm

**CODE** c-FOAM15



**POSITIONING  
FOAM,  
MATTRESS  
CUSTOM AS PER  
DRAWING**

Custom Mattress  
as per request

**CODE** c-  
FOAM16R



**POSITIONING  
FOAM, CUSTOM  
FOAM PIECE**

Custom Foam  
Piece: Ordered,  
as per request

**CODE** c-FOAM17



**POSITIONING  
FOAM BLOCK  
200X100X100M  
M**

Foam block with  
200x100x100mm

Square faces  
at each end.

**Foam Options:**

- Raw
- Covered in  
Dartex,
- Covered in  
Vinyl,
- Covered in  
Dartex  
Sealed.




**POSITIONING  
FOAM BLOCK  
180X180X600M  
M**

Foam positioning  
block with a  
180x180mm  
square face at  
each end.

**Foam Options:**

- Raw
  - Dartex
  - Vinyl
  - Sealed
- CODE** c-FOAM19




**POSITIONING  
FOAM SQUARE  
BLOCK  
250X250X50MM**

Foam Square  
block

**Foam Options:**

- Raw
  - Dartex
  - Vinyl
  - Sealed
- CODE** c-FOAM20







**POSITIONING  
FOAM TORSO  
BLOCK  
400X300X120M  
M**

Foam Torso block  
400x300x90mm  
deep.

**FOAM Options:**

- Raw
- or Covered in
- Dartex,
  - Vinyl,
- CODE** C-FOAM21

<div>CODE c-FOAM18</div>			
<div>  </div> <div> <b>POSITIONING SMALL ANIMAL CRADLE</b> </div> <div> Radiolucent positioning cradle for small, medium and large animals. </div> <div> <ul style="list-style-type: none"> <li>Small: 300mm long x 250mm x 120mm high</li> <li>Medium : 400mm long x 300mm x 120mm high</li> <li>Large: 500mm long x 350mm x 150mm high</li> </ul> </div> <div>CODE C-FOAM24</div>	<div>  </div> <div> <b>POSITIONING FOAM, INFANT IMMOBILIZER COVERED WITH DARTEX</b> </div> <div> Infant Immobilizer covered in DARTEX </div> <div> 1 x foam base (450x400x50 mm) with 2 x adjustable bolsters (450x100x90 mm) </div> <div> fastened and adjusted by velcro and covered with Dartex </div> <div>CODE c-FOAM60C</div>	<div>  </div> <div> <b>POSITIONING FOAM, HEAD IMMOBILISER 250X270X180M M</b> </div> <div> Head Immobiliser Foam: </div> <div> is a U shaped foam designed to support the head of a patient during X-ray, CT and MRI studies. 250x270x180mm </div> <div>CODE c-FOAM61</div>	<div>  </div> <div> <b>POSITIONING FOAM, FOREHEAD IMMOBILISER 120X270X250M M</b> </div> <div> Forehead Immobiliser: </div> <div> Is a U shaped designed to support the forehead of a patient during X-ray, CT and MRI studies. Can be covered in Dartex or vinyl. 120x270x250mm. </div> <div>CODE C-FOAM62</div>



**NECK INJECTION  
POSITIONER CT**

Nylon Neck  
Injection  
Positioning

Adjustable Strap.

For assisting patients to keep thier shoulders down during CT contrast injections. Strap is positioned under feet and two loops held tight to bring the shoulders down.

**CODE C-  
PRO1038**



**1KG SMALL  
SANDBAG**

**Double bagged  
Vinyl sandbags  
with a heavy-  
duty loop  
handle at one  
end for easy  
use.**

**CODE C-SA1001**



**3KG MEDIUM  
SANDBAG**

**Double bagged  
Vinyl sandbags  
with a heavy-  
duty loop  
handle at one  
end for easy  
use.**

**CODE C-SA1003**



**5KG LARGE  
SANDBAG**

**Double bagged  
Vinyl sandbags  
with a heavy-  
duty loop  
handle at one  
end for easy  
use.**

**CODE C-SA1005**



**SNAKE SANDBAG  
1.3KG**

100cm long x  
6cm wide - 1.3kg

Double bagged  
sandbag with a  
hard wearing  
Vinyl cover.

**CODE C-SA1006**



**SNAKE SANDBAG  
2.6KG**

100cm long x  
12cm wide -  
2.6kg

Double bagged  
sandbag with a  
hard wearing  
Vinyl cover.

**CODE C-SA1007**



**SANDBAG  
BUNDLE  
CONSISTING OF  
1KG, 3KG AND  
5KG**

**Sandbags  
double bagged  
with a loop  
handle at one  
end for easy  
use.**

**CODE C-SA1008**



**CUSTOM  
SANDBAG**

*Custom Vinyl  
Sandbag  
(double sealed)*

**CODE C-SA1009**



**RADIOGRAPHY  
CALIPER**

Easy to use  
caliper for  
anatomy  
thickness  
measurements.

Set you cm  
thickness  
generator factors  
for optimum  
plate dose.

**CODE E-  
ACC1033**



**CUSTOM EASY  
SLIDE, QUOTE  
ON REQUEST**

Custom Patient  
Easy Slide (non  
friction

material for easy  
sliding - made to  
order)

**CODE E-  
ACC1035**



**PATIENT EASY  
SLIDE  
LANDSCAPE**

Landscape Patient  
Easy Slide (non  
friction

material for easy  
sliding).

**CODE E-  
ACC1036**



**PATIENT EASY  
SLIDE PORTRAIT**

Portrait Patient  
Easy Slide (non  
friction material for  
easy sliding). Sold  
Individually.

**CODE E-  
ACC1037**