

COLLIMATOR FILTERS



**NOLAN FILTERS,  
A-P CERVICO  
THORACIC**

A-P Cervico  
Thoracic Compen  
sating wedge

**CODE** C-FIL1001



**NOLAN FILTER,  
A-P FEMALE  
GONAD SHIELD**

A-P Female  
Gonad Shield for  
AP shots.

**CODE** C-FIL1002



**NOLAN FILTER,  
A-P MALE  
GONAD**

A-P Male Gonad  
Shield for AP  
shots.

**CODE** C-FIL1003



**NOLAN FILTER,  
A-P PARASPINAL  
BLOCKER  
(TOTAL)**

A-P Paraspinal  
Blocker (Total):  
Totally blocks  
anatomy either  
side of the spine.

**CODE** C-FIL1004



**NOLAN FILTER,  
A-P PARASPINAL  
BLOCKER  
(PARTIAL)**

A-P Paraspinal  
Blocker (Partial):  
Partially blocks  
anatomy either  
side of the spine.


**CODE** C-FIL1005



**NOLAN FILTER,  
GONAD SHIELD  
SYSTEM SET**

Set of male and  
female gonad  
shields and belt.

**CODE** C-FIL1006



**NOLAN FILTER  
FACE PLATE**

Mounts directly  
on to the  
collimator to

accommodate  
multiple filters

**CODE** C-FIL1008



**NOLAN FILTER,  
LATERAL  
SPINOUS**

Lateral Spinous:  
Reduces scatter  
and skin dose  
and

balances plate  
dose by blocking  
x-rays to the  
spinous process.

**CODE** C-FIL1009



**NOLAN FILTER,  
GONAD SYSTEM  
BELT**

Velcro belt to  
suspend gonad  
shields

**CODE** C-FIL1010



**NOLAN FILTER  
LATERAL  
CERVICO  
THORACIC**

Accommodates  
relative thickness  
of the shoulder to  
the Cervical or  
Thyroid region.

**CODE** C-FIL1011



**NOLAN FILTER,  
LATERAL GONAD**

Lateral M/F  
Gonad shield for  
lateral exposures

**CODE** C-FIL1012



**NOLAN FILTER,  
LATERAL  
LUMBAR LUNG**

Lateral Lumbar  
Lung: Lumbar  
Lung shield for  
lateral exposures

**CODE** C-FIL1013



**NOLAN FILTER,  
LATERAL  
THORACIC**

Lateral Thoracic.  
Compensating  
wedge

**CODE** C-FIL1014



**NOLAN FILTER  
POINT FILTERS**

**POINT FILTERS.**  
*Can be used for  
compensation  
in all views.*

Adjust the  
exposure to the  
place according  
to anatomical  
thickness. Lower  
until the bottom  
is just above the  
iliac crests.

**CODE** C-FIL1015



**LATERAL  
STERNUM**

Lateral Sternum:  
Shield for lateral  
exposures

**CODE** C-FIL1017



**OCTOSTOP  
BOOMERANG  
FILTER**


Regular  
Boomerang L29  
x W19 x H4.5 cm  
approx.; 1.2 Kg.

**CODE** C-FIL1018



**OCTOSTOP  
INGOT FILTER**

Ingot Filter  
L29 x W19 x H4.5  
cm approx 1.8kg.  
**CODE** C-FIL1019



**OCTOSTOP HALF  
INGOT FILTER**

Half Ingot Filter  
L29 x W19 x H4.5  
cm approx 1.8kg.  
**CODE** C-FIL1034