

## COLLIMATOR FILTERS



**NOLAN FILTERS,  
A-P CERVICO  
THORACIC**

A-P Cervico  
Thoracic Compen  
sating wedge

**CODE NFAPCT**



**NOLAN FILTER,  
A-P FEMALE  
GONAD SHIELD**

A-P Female  
Gonad Shield for  
AP shots.

**CODE NFAPFGS**



**NOLAN FILTER,  
A-P MALE  
GONAD**

A-P Male Gonad  
Shield for AP  
shots.

**CODE NFAPMGS**



**NOLAN FILTER,  
GONAD SYSTEM  
BELT**

Velcro belt to  
suspend gonad  
shields

**CODE  
NFAPSBELT**



**NOLAN FILTER,  
A-P PARASPINAL  
BLOCKER  
(TOTAL)**

A-P Paraspinal  
Blocker (Total):  
Totally blocks  
anatomy either  
side of the spine.


**CODE NFAPSBF**



**NOLAN FILTER,  
A-P PARASPINAL  
BLOCKER  
(PARTIAL)**

A-P Paraspinal  
Blocker (Partial):  
Partially blocks  
anatomy either  
side of the spine.

**CODE NFAPSBP**




**NOLAN FILTER,  
GONAD SHIELD  
SYSTEM SET**

Set of male and  
female gonad  
shields

and belt.

**CODE NFAPSET**




**NOLAN FILTER  
FACE PLATE**

Mounts directly  
on to the  
collimator to

accommodate  
multiple filters

**CODE NFFP**



**FULL SET SPINE FILTERS**

Complete box set of all Full Spine Filters. **Reduce Radiation dose by up to 95%.**

**CODE NFFS**



**NOLAN FILTER LATERAL CERVICO THORACIC**

Accommodates relative thickness of the shoulder to the Cervical or Thyroid region.

**CODE NFLCT**



**NOLAN FILTER, LATERAL GONAD**

Lateral M/F Gonad shield for lateral exposures

**CODE NFLG**



**NOLAN FILTER, LATERAL LUMBAR LUNG**

Lateral Lumbar Lung: Lumbar Lung shield for lateral exposures

**CODE NFLLL**



**NOLAN FILTER, LATERAL SPINOUS**

Lateral Spinous: Reduces scatter and skin dose and

balances plate dose by blocking x-rays to the spinous process.

**CODE NFLS**



**NOLAN FILTER, LATERAL THORACIC**

Lateral Thoracic. Compensating wedge

**CODE NFLT**



**NOLAN FILTER POINT FILTERS**

**POINT FILTERS. Can be used for compensation in all views.**

Adjust the exposure to the plate according to anatomical thickness. Lower until the bottom is just above the iliac crests.

**CODE NFPF**



Lateral Sternum: Shield for lateral exposures

**CODE NFSB**



**OCTOSTOP  
BOOMERANG  
FILTER**

Regular  
Boomerang L29  
x W19 x H4.5 cm  
approx.; 1.2 Kg.  
**PRICE TO BE  
ADVICE**


**CODE** OCT01



**OCTOSTOP  
CHAINLESS  
BOOMERANG**

Regular  
Boomerang with  
no chain. **PRICE:  
TBA**

**CODE** OCT01NC




**OCTOSTOP  
INGOT FILTER**

Ingot Filter

L29 x W19 x H4.5  
cm approx 1.8kg.  
**PRICE: TBA**

**CODE** OCT02




**OCTOSTOP HALF  
INGOT FILTER**

Half Ingot Filter

L29 x W19 x H4.5  
cm approx 1.8kg.  
**PRICE: TBA**

**CODE** OCT02S




**OCTOSTOP  
GENTLE SLOPE  
FILTER**

Octostop Gentle  
Slope Filter.

L30x W17 x H2.3  
cm approx.; 800  
g. **PRICE: TBA**

**CODE** OCT03



**OCTOSTOP  
LITTLE PRISM  
FILTER**

Octostop Little  
Prism Filter.  
**PRICE: TBA**

**CODE** OCT04



**OCTOSTOP  
BOOMERANG  
FILTER 15 °**

Octostop  
Boomerang  
Filter 15 ° L21 x  
W7 x H3.3 cm.  
**PRICE: TBA**

**CODE** OCT06