

PERSONAL RADIATION PROTECTION



LEAD GOWN  
PRICE  
CALCULATOR

Choose your Size  
(XS) X Small: 38" / 96.5 cm Wide - Fits Chest 30-34" / 76.2-86.4 cm  
(S) Small: 20" / 51 cm Wide - Fits Chest 34-38" / 86.4-96.5 cm  
(M) Medium: 22" / 56 cm Wide - Fits Chest 38-42" / 96.5-106.7 cm  
(L) Large: 24" / 61 cm Wide - Fits Chest 42-46" / 106.7-116.8 cm

Regular lead =  
0.5mm

CODE C-  
PRO1000



OPEN PALM  
SHIELD

CODE C-  
PRO1001



VET MITT

CODE C-  
PRO1002



LEAD THYROID  
COLLAR

**Quick Release Collar** offers the same look and comfort as our model 621U, but without the hook & loop closure. Rear buckle closure. Available in our wide selection of nylons, ripstops and designer fabrics.

- **Standard Protection:** 0.5mm Pb Equivalency

**New Style Thyroid Collar** offers a trim yet shielded and comfortable design. Includes binding and rear hook & loop closure.

- **Standard Protection:** 0.5mm Pb Equivalency

Please note:  
Colour is light blue

CODE C-  
PRO1003-L



**WALL MOUNT  
SWING ARM  
APRON RACK**

Wall Mount Swing Arm Apron Rack features 5 individual chrome shoulder shape swing arms that fold flush against the wall. Available in right and left folding models. Simply hang two models side-by-side for a 10 hanger wall console unit.


**CODE C-  
PRO1005**



**PEG RACK WITH  
8 PEGS**

- Peg Rack is economical for the quick and easy storage of single or multiple aprons.
- Powder coated steel construction.
- Easy installation.

**CODE C-  
PRO1008**



**MOBILE APRON  
RACK 1.17MT**

**CODE C-  
PRO1009**



**HEAVY-DUTY  
APRON HANGER**

**CODE C-  
PRO1013**



**CLOSED LOOP  
APRON HANGER**

**CODE C-  
PRO1014**



**VINYL LEAD  
SEAMLESS  
GLOVES**

**CODE C-  
PRO1015**