

X-RAY COMPONENTS

DEL MEDICAL



DEL MEDICAL OVERHEAD TUBE CRANE SYSTEMS

Brochures and
Specifications

[OTC18T](#)

[OTC18M](#)

[VS300 Wall Stand](#)

[Table EV800
Specification](#)

CODE DMOTC



APIXIA EXL PSP SCANNER-108L

The EXL® PSP Scanner packs tremendous imaging power into a small package.

Tiny enough to fit into the most crowded work stations, yet robust enough to be placed in a central location for use by the entire clinic, the EXL® PSP Scanner is the ideal product for any size practice.

Included:

- 4 imaging plate for Size #2
- 1 imaging plate for Size #4
- 1 imaging plate for Size #5
- 300 pcs of plate barrier envelopes for Size #2
- 50 pcs of plate barrier envelopes for Size #4
- 50 pcs of plate barrier envelopes for Size #5
- Cleaning plates and Plate container
- USB cable and Ethernet cable
- Power adaptor, Imaging software and User manual

[View user
interface
demonstration
here](#)

[View Apixia EXL
PSP Scanner
Brochure here](#)

CODE E-
GEN1004



APIXIA EXM PSP SCANNER - 108M

The Apixia EXM® PSP Scanner packs tremendous imaging power into a small package. It is tiny enough to fit into the most crowded dental operatories for handy chairside use.

Imaging with Phosphor plates offers all the flexibility and patient comfort of film, combined with the immediacy and environmental benefits of digital radiography.

Included:

- 500 pcs of plate barrier envelopes for Size #2
- 12 imaging plates for Size #2 (3 boxes, one box has 4 pcs same size)
- Cleaning plates
- Plate container
- USB cable and Ethernet cable
- Power adaptor
- Imaging software and User manual
- Power adaptor, Imaging software and User manual

[View user
interface
demonstration
here](#)

[View Apixia ESM
PSP Scanner
Brochure here](#)

CODE E-
GEN1005



CPI CMP200 X- RAY GENERATOR

The **CMP 200® generator** series offers advanced high quality radiographic capabilities with unmatched performance and reliability resulting in the lowest cost of ownership over the product lifetime. This generator is suitable for film and CR-based radiographic systems and features leading edge computer-based controls to minimize patient dose and maximize image quality while demonstrating excellent reproducibility with user-friendly operator controls.

ATX has factory trained personnel to advise, install and maintain CPI generators.

Generators can be retrofitted into your existing bucky room or included into a new facility.

[CM Series Gen](#)

OPTIONS

- Automatic Exposure Control (AEC)
- Dose Area Product (DAP)
- Exposure Hand Switch
- Membrane Control Console
- Mini-Control Console

Lead time:

Subject to stock availability, transit times and radiation licensing. Estimated shipping date to be provided once the purchase is confirmed

CODE E-
GEN1016



MANUAL
COLLIMATOR
CML 125 -
COLLIMARE

Your First Choice
for Collimator
Needs

- Manual
Certified
Units
- Non-
Certified
Veterinari
an Units
- Portable
Certified
Units
- Portable
Non-
Certified
Veterinari
an Units
- Mobile
Certified
Units
- High
Throughp
ut
Clinical
Applicati
on Units
- Custom
Units
- Collimato
r
Accessori
es and
Add-ons



DEL FLOOR
MOUNTED TUBE
STAND

The **FMT Floor mounted tube stand** features an advanced, 10.4" technologist-friendly touchscreen console that provides large text readouts and integrated operation of the tube stand. Ergonomically placed movement identification buttons provide optimal accuracy when setting up views. Take advantage of extensive longitudinal and vertical movements, as well as a 180° rotating tube column enabling positioning for cross-table lateral exams.

The FMT Floor Mounted Tube Stand comes in 3 different models:

[FMT18T](#)

[FMT18M](#)

[FMT](#)

CODE E-
GEN1032



DEL RT100
PATIENT TABLE

The sturdy **RT100 four-way float top table** is suitable for use in a variety of general radiographic environments including small hospitals, clinics, ambulatory care centers and private offices. The table features a 700 lbs (318 kg) patient weight capacity and a stain resistant floating table top which glides freely on a precision roller bearing system. Float top movements are safely controlled by heavy-duty electric locks, activated by a conveniently placed foot treadle release. The 17" x 17" standard bucky tray can accommodate any cassette sized detector with ease. Its all steel structure and welded base provides maximum strength and reliability. In addition, the RT100's pedestal base can be used to store generator electronics making it a great space



X-RAY TUBE
VAREX
DIAMOND RAD
14

The RAD-14 is a 3" (80mm), 150 kV, 212 kJ (300 KHU) rotating anode insert specifically designed for general radiographic and fluoro/spotfilm procedures. The insert features a 12° tungsten rhenium molybdenum target and is available in the following focal spot combinations:

0.3 - 1.2
0.6 - 1.2
0.6 - 1.5
IEC 60336

Nominal Anode
Input Power:
Small (0.3) - 7.5
kW IEC 60613
Small (0.6) - 32
kW IEC 60613
Large (1.2) - 77
kW IEC 60613
Large (1.5) - 95
kW IEC 60613For
the equivalent
anode input
power of 80
Watts

Lead time:
Subject to stock availability, transit times and radiation licensing. Estimated shipping date to be provided once



Lead time:
Subject to stock availability, transit times and radiation licensing. Estimated shipping date to be provided once the purchase is confirmed

CODE E-GEN1020

saver.

[RT100](#)

[Please contact us for pricing](#)

CODE E-GEN1041

the purchase is confirmed

CODE E-GEN1083



IRAY CANIS DENTAL CAMERA X RAY

The IRay Portable Camera XRay is a portable intraoral X-ray device.

It is equipped with a 7-inch LCD touch screen and Android operating system. It can utilise an iRay sensor, and images can be transferred directly via WiFi or USB, allowing for simple point-and-shoot functionality. The machine is compact and lightweight, with a variety of network connection options.

CODE E-GEN1117



ATX PORTABLE X-RAY GENERATOR 5KW

ATX Portable X-ray system


- Designed with innovation in mind, Ultra-light, high frequency x-ray system.
- Compact size and lightweight body
- Durable rugged exterior

Provided with Compliance testing as per local radiation health regulations.

Lead time:
Subject to stock availability, transit times and radiation licensing. Estimated shipping date to be provided once the purchase is confirmed

[Please contact us for pricing](#)

CODE E-GEN1123



ATX PORTABLE X-RAY GENERATOR 2.4KW

ATX Portable X-ray system

- Designed with innovation in mind, Ultra-light, high frequency x-ray system.
- Compact size and lightweight body
- Durable rugged exterior

Provided with Compliance testing as per local radiation health regulations.

Lead time:
Subject to stock availability, transit times and radiation licensing. Estimated shipping date to be provided once the purchase is confirmed

[Please contact us for pricing](#)

CODE E-GEN1124



ATX PORTABLE X-RAY GENERATOR 1.8KW

ATX Portable X-ray system

- Designed with innovation in mind, Ultra-light, high frequency x-ray system.
- Compact size and lightweight body
- Durable rugged exterior

Provided with Compliance testing as per local radiation health regulations.

Lead time:
Subject to stock availability, transit times and radiation licensing. Estimated shipping date to be provided once the purchase is confirmed

[Please contact us for pricing](#)

CODE E-GEN1125



ATX SMART
PORTABLE X-
RAY GENERATOR

ATX Smart-
portable X-ray
system

- Embedde
d tablet
touch
screen
console
- Designed
with
innovatio
n in
mind,
Ultra-
light,
high
frequenc
y x-ray
system.
- Compact
size and
lightweig
ht body
- Durable
rugged
exterior
- Auto
calibratio
n
function
- Body part
APR
function
(over 960
APR)
- DAP
display
function
- Also
available
with a
traditiona
l button
keypad
control

Compliance
testing as per
local radiation
health
regulations
proivided

Lead time:
Subject to stock
availability,
transit times and
radiation
licensing.
Estimated
shipping date to
be provided once
the purchase is
confirmed

CODE EQU-
GEN-10020



PODOBLOCK. 90
CODE
PPB01AL090



PODOBLOCK.
120
CODE
PPB01AL120



PODOBLOCK.
170
CODE
PPB01AL170

