

## X-RAY COMPONENTS

### DEL MEDICAL



#### DEL MEDICAL OVERHEAD TUBE CRANE SYSTEMS

Brochures and  
Specifications

[OTC18T](#)

[OTC18M](#)

[VS300 Wall Stand](#)

[Table EV800  
Specification](#)

**CODE** DMOTC



#### APIXIA EXL PSP SCANNER-108L

The EXL® PSP Scanner packs tremendous imaging power into a small package.

Tiny enough to fit into the most crowded work stations, yet robust enough to be placed in a central location for use by the entire clinic, the EXL® PSP Scanner is the ideal product for any size practice.

##### **Included:**

- 4 imaging plate for Size #2
- 1 imaging plate for Size #4
- 1 imaging plate for Size #5
- 300 pcs of plate barrier envelopes for Size #2
- 50 pcs of plate barrier envelopes for Size #4
- 50 pcs of plate barrier envelopes for Size #5
- Cleaning plates and Plate container
- USB cable and Ethernet cable
- Power adaptor, Imaging software and User manual

[View user interface demonstration here](#)

[View Apixia EXL PSP Scanner Brochure here](#)

**CODE** E-  
GEN1004



#### APIXIA EXM PSP SCANNER - 108M

**The Apixia EXM® PSP Scanner** packs tremendous imaging power into a small package. It is tiny enough to fit into the most crowded dental operatories for handy chairside use.

Imaging with Phosphor plates offers all the flexibility and patient comfort of film, combined with the immediacy and environmental benefits of digital radiography.

##### **Included:**

- 500 pcs of plate barrier envelopes for Size #2
- 12 imaging plates for Size #2 (3 boxes, one box has 4 pcs same size)
- Cleaning plates
- Plate container
- USB cable and Ethernet cable
- Power adaptor
- Imaging software and User manual
- Power adaptor, Imaging software and User manual

[View user interface demonstration here](#)

[View Apixia ESM PSP Scanner Brochure here](#)

**CODE** E-  
GEN1005



#### CPI CMP200 X- RAY GENERATOR

The **CMP 200® generator** series offers advanced high quality radiographic capabilities with unmatched performance and reliability resulting in the lowest cost of ownership over the product lifetime. This generator is suitable for film and CR-based radiographic systems and features leading edge computer-based controls to minimize patient dose and maximize image quality while demonstrating excellent reproducibility with user-friendly operator controls.

**ATX has factory trained personnel to advise, install and maintain CPI generators.**

Generators can be retrofitted into your existing bucky room or included into a new facility.

[CM Series Gen](#)

##### **OPTIONS**

- Automatic Exposure Control (AEC)
- Dose Area Product (DAP)
- Exposure Hand Switch
- Membrane Control Console
- Mini-Control Console

**Lead time:**

Subject to stock availability, transit times and radiation licensing. Estimated shipping date to be provided once the purchase is confirmed

CODE E-  
GEN1016



MANUAL  
COLLIMATOR  
CML 125 -  
COLLIMARE

Your First Choice  
for Collimator  
Needs

- Manual  
Certified  
Units
- Non-  
Certified  
Veterinari  
an Units
- Portable  
Certified  
Units
- Portable  
Non-  
Certified  
Veterinari  
an Units
- Mobile  
Certified  
Units
- High  
Throughp  
ut  
Clinical  
Applicati  
on Units
- Custom  
Units
- Collimato  
r  
Accessori  
es and  
Add-ons



DEL FLOOR  
MOUNTED TUBE  
STAND

The **FMT Floor mounted tube stand** features an advanced, 10.4" technologist-friendly touchscreen console that provides large text readouts and integrated operation of the tube stand. Ergonomically placed movement identification buttons provide optimal accuracy when setting up views. Take advantage of extensive longitudinal and vertical movements, as well as a 180° rotating tube column enabling positioning for cross-table lateral exams.

The FMT Floor Mounted Tube Stand comes in 3 different models:

[FMT18T](#)

[FMT18M](#)

[FMT](#)

CODE E-  
GEN1032



DEL RT100  
PATIENT TABLE

The sturdy **RT100 four-way float top table** is suitable for use in a variety of general radiographic environments including small hospitals, clinics, ambulatory care centers and private offices. The table features a 700 lbs (318 kg) patient weight capacity and a stain resistant floating table top which glides freely on a precision roller bearing system. Float top movements are safely controlled by heavy-duty electric locks, activated by a conveniently placed foot treadle release. The 17" x 17" standard bucky tray can accommodate any cassette sized detector with ease. Its all steel structure and welded base provides maximum strength and reliability. In addition, the RT100's pedestal base can be used to store generator electronics making it a great space



X-RAY TUBE  
VAREX  
DIAMOND RAD  
14

The RAD-14 is a 3" (80mm), 150 kV, 212 kJ (300 KHU) rotating anode insert specifically designed for general radiographic and fluoro/spotfilm procedures. The insert features a 12° tungsten rhenium molybdenum target and is available in the following focal spot combinations:

0.3 - 1.2  
0.6 - 1.2  
0.6 - 1.5  
IEC 60336

Nominal Anode  
Input Power:  
Small (0.3) - 7.5  
kW IEC 60613  
Small (0.6) - 32  
kW IEC 60613  
Large (1.2) - 77  
kW IEC 60613  
Large (1.5) - 95  
kW IEC 60613For  
the equivalent  
anode input  
power of 80  
Watts

**Lead time:**  
Subject to stock availability, transit times and radiation licensing. Estimated shipping date to be provided once



**Lead time:**  
Subject to stock availability, transit times and radiation licensing. Estimated shipping date to be provided once the purchase is confirmed

**CODE E-GEN1020**

saver.

[RT100](#)

[Please contact us for pricing](#)

**CODE E-GEN1041**

the purchase is confirmed

**CODE E-GEN1083**



**IRAY CANIS DENTAL CAMERA X RAY**

The IRay Portable Camera XRay is a portable intraoral X-ray device.

It is equipped with a 7-inch LCD touch screen and Android operating system. It can utilise an iRay sensor, and images can be transferred directly via WiFi or USB, allowing for simple point-and-shoot functionality. The machine is compact and lightweight, with a variety of network connection options.

**CODE E-GEN1117**



**ATX PORTABLE X-RAY GENERATOR 5KW**

**ATX Portable X-ray system**


- Designed with innovation in mind, Ultra-light, high frequency x-ray system.
- Compact size and lightweight body
- Durable rugged exterior

Provided with Compliance testing as per local radiation health regulations.

**Lead time:**  
Subject to stock availability, transit times and radiation licensing. Estimated shipping date to be provided once the purchase is confirmed

[Please contact us for pricing](#)

**CODE E-GEN1123**



**ATX PORTABLE X-RAY GENERATOR 2.4KW**

**ATX Portable X-ray system**

- Designed with innovation in mind, Ultra-light, high frequency x-ray system.
- Compact size and lightweight body
- Durable rugged exterior

Provided with Compliance testing as per local radiation health regulations.

**Lead time:**  
Subject to stock availability, transit times and radiation licensing. Estimated shipping date to be provided once the purchase is confirmed

[Please contact us for pricing](#)

**CODE E-GEN1124**



**ATX PORTABLE X-RAY GENERATOR 1.8KW**

**ATX Portable X-ray system**

- Designed with innovation in mind, Ultra-light, high frequency x-ray system.
- Compact size and lightweight body
- Durable rugged exterior

Provided with Compliance testing as per local radiation health regulations.

**Lead time:**  
Subject to stock availability, transit times and radiation licensing. Estimated shipping date to be provided once the purchase is confirmed

[Please contact us for pricing](#)

**CODE E-GEN1125**



ATX SMART  
PORTABLE X-  
RAY GENERATOR

ATX Smart-  
portable X-ray  
system

- Embedde  
d tablet  
touch  
screen  
console
- Designed  
with  
innovatio  
n in  
mind,  
Ultra-  
light,  
high  
frequenc  
y x-ray  
system.
- Compact  
size and  
lightweig  
ht body
- Durable  
rugged  
exterior
- Auto  
calibratio  
n  
function
- Body part  
APR  
function  
(over 960  
APR)
- DAP  
display  
function
- Also  
available  
with a  
traditiona  
l button  
keypad  
control

Compliance  
testing as per  
local radiation  
health  
regulations  
proivided

**Lead time:**  
Subject to stock  
availability,  
transit times and  
radiation  
licensing.  
Estimated  
shipping date to  
be provided once  
the purchase is  
confirmed

**CODE** EQU-  
GEN-10020



**PODOBLOCK. 90**  
**CODE**  
PPB01AL090



**PODOBLOCK. 120**  
**CODE**  
PPB01AL120



**PODOBLOCK. 170**  
**CODE**  
PPB01AL170

