

# X-RAY EQUIPMENT



## ATX ULTRA 99V DR PACKAGE

**Highest resolution, 99 micron, accurate imaging.**

The ULTRA-99V is a cassette-sized, next generation flat panel DR detector. Battery stands by for atleast 12 hours and supports IEEE 802.11ac with transmission speed up to 1300Mbps. The detector acquires outstanding images with high DQE and high resolutions. Equipped with highly sensitive AED function, the detector can be easily connected to and synchronized with any kind of X-ray generator

**Included on package:**

- ATX Flat DR Panel - 14x17inch - Wireless Csl scintillator
- Lenovo ThinkCentre M920z All-in-One Computer - Core i5 i5-8500 - 8 GB RAM - 512 GB SSD - 60.5 cm (23.8") 1920 x 1080 - Desktop
- VET DROC Acquisition Software
- Training & Installation
- iTX Cloud Pacs



ATX DR  
ULTRA-140H

## ULTRA-140H DR PANEL PACKAGE

**Included in DR Package:**

- ATX Flat DR Panel - 10x12inch  
14x17inch  
17x17inch - Wireless or Wired / Csl or Gadox
- Lenovo Touch screen all in one computer
- Ecom DROC Acquisition Software
- Training & Installation
- ATX Cloud Pacs DICOM Image Storage platform
- 2 Years warranty & Assurance support

**CODE EQU-GEN-10010**



DICOM  
Image  
Storage  
platform-  
one year  
free trial

- 2 Years  
standard  
warranty

**Lead time:**

Subject to stock  
availability, transit  
times and radiation  
licensing.  
Estimated shipping  
date to be provided  
once the purchase  
is confirmed

**CODE** E-  
GEN1102

