

MARKERS



**AKUCAL
SCALING TOOL (25MM
REPLACEMENT
MARKER)**

25mm

Replacement
stainless steel
marker for Image
Scaling Tool

EACH

**CODE E-
ACC1030-25**




**AKUCAL IMAGE
SCALING DEVICE
(30MM)
(COMPLETE)**

25mm image
scaling device

COMPLETE with
all parts
including-

- Suction
cup base for
Horizontal or
Vertical
use
- Flexible
arm to
position
marker at
anatomy
level
- Aluminium
tip
avoids
breakage
from
wear and
tear
- Numerical
size
indicators
visible in
every
image (25mm,
25.4mm,
30mm)
- Surgical
Grade
stainless
steel
marker
- Lever to
engage
or
disengage
suction

**CODE E-
ACC1031-30**



**AKUCAL
SCALING TOOL (30MM
REPLACEMENT
MARKER)**

30mm

Replacement
stainless steel
marker for Image
Scaling Tool

EACH

**CODE E-
ACC1032-30**

