

X-RAY EQUIPMENT



ATX ULTRA 99V
DR PACKAGE

Highest resolution, 99 micron, accurate imaging.

The ULTRA-99V is a cassette-sized, next generation flat panel DR detector. Battery stands by for atleast 12 hours and supports IEEE 802.11ac with transmission speed up to 1300Mbps. The detector acquires outstanding images with high DQE and high resolutions. Equipped with highly sensitive AED function, the detector can be easily connected to and synchronized with any kind of X-ray generator

Included on package:

- ATX Flat DR Panel - 14x17inch - Wireless Csl scintillator
- Lenovo ThinkCentre M920z All-in-One Computer - Core i5 i5-8500 - 8 GB RAM - 512 GB SSD - 60.5 cm (23.8") 1920 x 1080 - Desktop
- VET DROC Acquisition Software
- Training & Installation
- iTX Cloud Pacs



ATX DR
ULTRA-140H

ULTRA-140H DR
PANEL PACKAGE

Included in DR Package:

- ATX Flat DR Panel - 10x12inch 14x17inch 17x17inch - Wireless or Wired / Csl or Gadox
- Lenovo Touch screen all in one computer
- Ecom DROC Acquisition Software
- Training & Installation
- ATX Cloud Pacs DICOM Image Storage platform
- 2 Years warranty & Assurance support

CODE EQU-GEN-10010



DICOM
Image
Storage
platform-
one year
free trial

- 2 Years
standard
warranty

Lead time:

Subject to stock
availability, transit
times and radiation
licensing.
Estimated shipping
date to be provided
once the purchase
is confirmed

CODE E-
GEN1102

