

MARKERS




**AKUCAL IMAGE  
SCALING DEVICE  
( 25MM )  
(COMPLETE)**

25mm image  
scaling device

COMPLETE with  
all parts  
including-

- Suction cup base for Horizontal or Vertical use
- Flexible arm to position marker at anatomy level
- Aluminium tip avoids breakage from wear and tear
- Numerical size indicators visible in every image ( 25mm, 25.4mm, 30mm )
- Surgical Grade stainless steel marker
- Lever to engage or disengage suction

**CODE E-  
ACC1029-25**



**AKUCAL  
SCALING TOOL ( 25MM  
REPLACEMENT  
MARKER )**

25mm

Replacement  
stainless steel  
marker for Image  
Scaling Tool

EACH

**CODE E-  
ACC1030-25**



**AKUCAL  
SCALING TOOL ( 30MM  
REPLACEMENT  
MARKER )**

30mm

Replacement  
stainless steel  
marker for Image  
Scaling Tool

EACH

**CODE E-  
ACC1032-30**

