

CR AND DR ACCESSORIES



ATX CR PLATE CLEANER

ATX CR Plate Cleaner, is easy to use; Plate Cleaning Solution 250mls per bottle.

CR X-ray Plates are highly sensitive and requires regular cleaning. We recommend cleaning your plate every 200 scans, or every 30 days (depending of your work environment).

CODE SC300ATX



ATX CR Screen Cleaner Kit

ATX CR PLATE CLEANER KIT

ATX CR Plate Cleaner Kit:

- 1 x 250ml bottle ATX CR Screen Cleaner
- 1 x Air Bulb Brush
- 1 x ATX Equipment Wipes (75 pack)

CODE SC300KIT

