

MARKERS



BEEKLEY X-SPOT® 101

BEEKLEY X-SPOT 101 BLUE

Ref: 101 Blue,
(QTY: 150/box)

X-SPOT 101: Identify and differentiate specific points of reference in conventional simulation treatment planning, chest, and trauma imaging.

- 1.5mm lead-free pellet
- easy Pick-Up tab

CODE C-
BEEK101



BEEKLEY Y-SPOT® 102

BEEKLEY Y-SPOT 102 Green

Ref: 102 Green
(QTY: 132/box)

Identify and differentiate specific points of reference in conventional simulation treatment planning, chest, and trauma imaging.

- 2.0mm lead-free pellet
- easy Pick-Up tab

CODE C-
BEEK102



BEEKLEY D-SPOT® ARROW 104

BEEKLEY D-SPOT® ARROW 104

Ref: 104 (QTY: 99/box)

- lead-free arrow with easy Pick-Up tab
- for marking nodes and outer canthus

CODE C-
BEEK104



BEEKLEY T-SPOT® LINE 111

BEEKLEY T-SPOT® LINE 111

Ref: 111, (QTY: 275cm cut to measure roll/box)

- .8mm diameter metallic line.
- for marking scars, nodes, field borders, and match lines

CODE C-
BEEK111

