

## POSITIONING AIDS



### TMJ POSITIONING BOX OF 24

TMJ (temporomandibular joint) positioning wedge box of 24 are invaluable for detailed insights into TMJ disorders. It helps in positioning the patient's head and jaw accurately for imaging procedures like X-rays or CT scans of the TMJ.

#### Box of 24 Clickers:

- Eliminates broken clickers during use
- User-friendly for patients
- Adjusts in 4mm increments with a total opening of 8 steps or 40mm
- Disposable after each use

These clickers can be cleaned like any other tool you use—through autoclaving, cleaning solutions, etc. They are shipped non-sterile and require cleaning before use.

**CODE** c-gen1060



