

CONTRASTS



BARIUM PASTE 20ML

Barium Sulphate, radiopaque paste that provides excellent radiographic visualization of the true hoof wall. The tube contains 20mls of paste. The paste appears as an opaque (white line) on the image, when applied in a continuous bead from the coronary band to the ground surface. Due to the very thin depth of horn at the most proximal crown of the hoof wall the actual horn margin is rarely visible on even the softest penetrating views.

**CODE C-
GEN1010**

MARKERS




BARIUM PASTE 20ML

Barium Sulphate, radiopaque paste that provides excellent radiographic visualization of the true hoof wall. The tube contains 20mls of paste. The paste appears as an opaque (white line) on the image, when applied in a continuous bead from the coronary band to the ground surface. Due to the very thin depth of horn at the most proximal crown of the hoof wall the actual horn margin is rarely visible on even the softest penetrating views.

**CODE C-
GEN1010**

X-RAY EQUIPMENT



**ATX PORTABLE
X-RAY
GENERATOR
5KW**

ATX Portable X-ray system

- Designed with innovation in mind, Ultra-light, high frequency x-ray system.
- Compact size and lightweight body
- Durable rugged exterior

Provided with Compliance testing as per local radiation health regulations.

Lead time:
Subject to stock availability, transit times and radiation licensing. Estimated shipping date to be provided once the purchase is confirmed

[Please contact us for pricing](#)

**CODE E-
GEN1123**



**ATX PORTABLE
X-RAY
GENERATOR
2.4KW**

ATX Portable X-ray system

- Designed with innovation in mind, Ultra-light, high frequency x-ray system.
- Compact size and lightweight body
- Durable rugged exterior

Provided with Compliance testing as per local radiation health regulations.

Lead time:
Subject to stock availability, transit times and radiation licensing. Estimated shipping date to be provided once the purchase is confirmed

[Please contact us for pricing](#)

**CODE E-
GEN1124**



**ATX PORTABLE
X-RAY
GENERATOR
1.8KW**

ATX Portable X-ray system

- Designed with innovation in mind, Ultra-light, high frequency x-ray system.
- Compact size and lightweight body
- Durable rugged exterior

Provided with Compliance testing as per local radiation health regulations.

Lead time:
Subject to stock availability, transit times and radiation licensing. Estimated shipping date to be provided once the purchase is confirmed

[Please contact us for pricing](#)

**CODE E-
GEN1125**

