

# COLLIMATOR FILTERS



**FULL SET SPINE FILTERS**

Complete box set of all Full Spine Filters. **Reduce Radiation dose by up to 95%.**

**CODE** C-FIL1025



**OCTOSTOP CHAINLESS BOOMERANG**

Regular Boomerang with no chain. **PRICE: TBA**

**CODE** C-FIL1033

**RADIATION PROTECTION**



**NOLAN FILTER,  
A-P FEMALE  
GONAD SHIELD**

A-P Female  
Gonad Shield for  
AP shots.

**CODE C-FIL1002**



**NOLAN FILTER,  
A-P MALE  
GONAD**

A-P Male Gonad  
Shield for AP  
shots.

**CODE C-FIL1003**



**NOLAN FILTER,  
A-P PARASPINAL  
BLOCKER  
(TOTAL)**

A-P Paraspinal  
Blocker (Total):  
Totally blocks  
anatomy either  
side of the spine.

**CODE C-FIL1004**



**NOLAN FILTER,  
A-P PARASPINAL  
BLOCKER  
(PARTIAL)**

A-P Paraspinal  
Blocker (Partial):  
Partially blocks  
anatomy either  
side of the spine.

**CODE C-FIL1005**



**NOLAN FILTER,  
LATERAL  
SPINOUS**

Lateral Spinous:  
Reduces scatter  
and skin dose  
and

balances plate  
dose by blocking  
x-rays to the  
spinous process.

**CODE C-FIL1009**



**NOLAN FILTER,  
GONAD SYSTEM  
BELT**

Velcro belt to  
suspend gonad  
shields

**CODE C-FIL1010**



**NOLAN FILTER  
LATERAL  
CERVICO  
THORACIC**

Accommodates  
relative thickness  
of the shoulder to  
the Cervical or  
Thyroid region.

**CODE C-FIL1011**



**NOLAN FILTER,  
LATERAL GONAD**

Lateral M/F  
Gonad shield for  
lateral exposures

**CODE C-FIL1012**



**NOLAN FILTER,  
LATERAL  
LUMBAR LUNG**

Lateral Lumbar Lung: Lumbar Lung shield for lateral exposures

**CODE C-FIL1013**



**NOLAN FILTER,  
LATERAL  
THORACIC**

Lateral Thoracic. Compensating wedge

**CODE C-FIL1014**



**LATERAL  
STERNUM**

Lateral Sternum: Shield for lateral exposures

**CODE C-FIL1017**



**OCTOSTOP  
BOOMERANG  
FILTER**

Regular Boomerang L29 x W19 x H4.5 cm approx.; 1.2 Kg.

**CODE C-FIL1018**




**OCTOSTOP  
INGOT FILTER**

Ingot Filter

L29 x W19 x H4.5 cm approx 1.8kg.

**CODE C-FIL1019**



**OCTOSTOP HALF  
INGOT FILTER**

Half Ingot Filter

L29 x W19 x H4.5 cm approx 1.8kg.

**CODE C-FIL1034**

