


COLLIMATOR FILTERS



**FULL SET SPINE
FILTERS**

Complete box
set of all Full
Spine
Filters. **Reduce
Radiation dose
by up to 95%.**

CODE C-FIL1025



**OCTOSTOP
CHAINLESS
BOOMERANG**

Regular
Boomerang with
no chain. **PRICE:
TBA**

CODE C-FIL1033

RADIATION PROTECTION



**NOLAN FILTER,
A-P FEMALE
GONAD SHIELD**

A-P Female
Gonad Shield for
AP shots.

CODE C-FIL1002



**NOLAN FILTER,
A-P MALE
GONAD**

A-P Male Gonad
Shield for AP
shots.

CODE C-FIL1003



**NOLAN FILTER,
A-P PARASPINAL
BLOCKER
(TOTAL)**

A-P Paraspinal
Blocker (Total):
Totally blocks
anatomy either
side of the spine.

CODE C-FIL1004



**NOLAN FILTER,
A-P PARASPINAL
BLOCKER
(PARTIAL)**

A-P Paraspinal
Blocker (Partial):
Partially blocks
anatomy either
side of the spine.

CODE C-FIL1005



**NOLAN FILTER,
LATERAL
SPINOUS**

Lateral Spinous:
Reduces scatter
and skin dose
and

balances plate
dose by blocking
x-rays to the
spinous process.

CODE C-FIL1009



**NOLAN FILTER,
GONAD SYSTEM
BELT**

Velcro belt to
suspend gonad
shields

CODE C-FIL1010



**NOLAN FILTER
LATERAL
CERVICO
THORACIC**

Accommodates
relative thickness
of the shoulder to
the Cervical or
Thyroid region.

CODE C-FIL1011



**NOLAN FILTER,
LATERAL GONAD**

Lateral M/F
Gonad shield for
lateral exposures

CODE C-FIL1012



**NOLAN FILTER,
LATERAL
LUMBAR LUNG**

Lateral Lumbar Lung: Lumbar Lung shield for lateral exposures

CODE C-FIL1013



**NOLAN FILTER,
LATERAL
THORACIC**

Lateral Thoracic: Compensating wedge

CODE C-FIL1014



**LATERAL
STERNUM**

Lateral Sternum: Shield for lateral exposures

CODE C-FIL1017



**OCTOSTOP
BOOMERANG
FILTER**

Regular Boomerang L29 x W19 x H4.5 cm approx.; 1.2 Kg.

CODE C-FIL1018




**OCTOSTOP
INGOT FILTER**

Ingot Filter

L29 x W19 x H4.5 cm approx 1.8kg.

CODE C-FIL1019



**OCTOSTOP HALF
INGOT FILTER**

Half Ingot Filter

L29 x W19 x H4.5 cm approx 1.8kg.

CODE C-FIL1034

