


COLLIMATOR FILTERS



**OCTOSTOP
BOOMERANG
FILTER 15 °**

Octostop
Boomerang
Filter 15 ° L21 x
W7 x H3.3 cm.
PRICE: TBA

CODE C-FIL1022



**FULL SET SPINE
FILTERS**

Complete box
set of all Full
Spine
Filters. **Reduce
Radiation dose
by up to 95%.**

CODE C-FIL1025

RADIATION PROTECTION



**NOLAN FILTERS,
A-P CERVICO
THORACIC**

A-P Cervico
Thoracic Compen
sating wedge

CODE C-FIL1001



**NOLAN FILTER,
A-P FEMALE
GONAD SHIELD**

A-P Female
Gonad Shield for
AP shots.

CODE C-FIL1002



**NOLAN FILTER,
A-P MALE
GONAD**

A-P Male Gonad
Shield for AP
shots.

CODE C-FIL1003



**NOLAN FILTER,
A-P PARASPINAL
BLOCKER
(TOTAL)**

A-P Paraspinal
Blocker (Total):
Totally blocks
anatomy either
side of the spine.


CODE C-FIL1004



**NOLAN FILTER,
A-P PARASPINAL
BLOCKER
(PARTIAL)**

A-P Paraspinal
Blocker (Partial):
Partially blocks
anatomy either
side of the spine.

CODE C-FIL1005




**NOLAN FILTER,
GONAD SHIELD
SYSTEM SET**

Set of male and
female gonad
shields

and belt.

CODE C-FIL1006



**NOLAN FILTER
FACE PLATE**

Mounts directly
on to the
collimator to

accommodate
multiple filters

CODE C-FIL1008



**NOLAN FILTER,
LATERAL
SPINOUS**

Lateral Spinous:
Reduces scatter
and skin dose
and

balances plate
dose by blocking
x-rays to the
spinous process.

CODE C-FIL1009



**NOLAN FILTER,
GONAD SYSTEM
BELT**

Velcro belt to
suspend gonad
shields

CODE C-FIL1010



**NOLAN FILTER
LATERAL
CERVICO
THORACIC**

Accommodates
relative thickness
of the shoulder to
the Cervical or
Thyroid region.

CODE C-FIL1011



**NOLAN FILTER,
LATERAL GONAD**

Lateral M/F
Gonad shield for
lateral exposures

CODE C-FIL1012



**NOLAN FILTER,
LATERAL
LUMBAR LUNG**

Lateral Lumbar
Lung: Lumbar
Lung shield for
lateral exposures

CODE C-FIL1013



**NOLAN FILTER,
LATERAL
THORACIC**

Lateral Thoracic.
Compensating
wedge

CODE C-FIL1014



**NOLAN FILTER
POINT FILTERS**

POINT FILTERS.
*Can be used for
compensation
in all views.*

Adjust the
exposure to the
place according
to anatomical
thickness. Lower
until the bottom
is just above the
iliac crests.

CODE C-FIL1015



**LATERAL
STERNUM**

Lateral Sternum:
Shield for lateral
exposures

CODE C-FIL1017

