

## X-RAY EQUIPMENT



### DEL FLOOR MOUNTED TUBE STAND

The **FMT Floor mounted tube stand** features an advanced, 10.4" technologist-friendly touchscreen console that provides large text readouts and integrated operation of the tube stand. Ergonomically placed movement identification buttons provide optimal accuracy when setting up views. Take advantage of extensive longitudinal and vertical movements, as well as a 180° rotating tube column enabling positioning for cross-table lateral exams.

The FMT Floor Mounted Tube Stand comes in 3 different models:

[FMT18T](#)

[FMT18M](#)

[FMT](#)

**CODE E-  
GEN1032**

