


COLLIMATOR FILTERS



FULL SET SPINE FILTERS

Complete box
set of all Full
Spine
Filters. **Reduce
Radiation dose
by up to 95%.**

CODE C-FIL1025


RADIATION PROTECTION



**NOLAN FILTERS,
A-P CERVICO
THORACIC**

A-P Cervico
Thoracic Compen
sating wedge

CODE C-FIL1001



**NOLAN FILTER
FACE PLATE**

Mounts directly
on to the
collimator to

accommodate
multiple filters

CODE C-FIL1008



**NOLAN FILTER
POINT FILTERS**

***POINT FILTERS.
Can be used for
compensation
in all views.***

Adjust the
exposure to the
place according
to anatomical
thickness. Lower
until the bottom
is just above the
iliac crests.

CODE C-FIL1015

