

## COLLIMATOR FILTERS



### FULL SET SPINE FILTERS

Complete box set of all Full Spine Filters. **Reduce Radiation dose by up to 95%.**

**CODE** C-FIL1025


# RADIATION PROTECTION



**NOLAN FILTERS,  
A-P CERVICO  
THORACIC**

A-P Cervico  
Thoracic Compen  
sating wedge

**CODE C-FIL1001**



**NOLAN FILTER  
FACE PLATE**

Mounts directly  
on to the  
collimator to

accommodate  
multiple filters

**CODE C-FIL1008**



**NOLAN FILTER  
POINT FILTERS**

***POINT FILTERS.  
Can be used for  
compensation  
in all views.***

Adjust the  
exposure to the  
place according  
to anatomical  
thickness. Lower  
until the bottom  
is just above the  
iliac crests.

**CODE C-FIL1015**

