

MARKERS



BEEKLEY
GUIDELINES® CT
BIOPSY GRID
117

GuideLines® CT
Biopsy
Grid 117

Ref: 117,
(Qty: 20/box)

- 4" x 5" biopsy grid
- 7 non-metallic radiopaque lines
- open pane design for easy mark placement
- ideal for CT-guided biopsies of the neck and other small areas

CODE C-
BEEK117



BEEKLEY
GUIDELINES® CT
BIOPSY GRID
217

GuideLines® CT
Biopsy
Grid 217

Ref: 217,
(Qty: 15/box)

- 7" x 7.5" biopsy grid
- 12 non-metallic radiopaque lines
- open pane design for easy mark placement
- ideal for CT-guided biopsies of abdomen and other large areas

CODE C-
BEEK217

