

INFECTION CONTROL



DISPOSABLE PATIENT GOWN (BOX OF 100)

A simple change can have a large impact on patient dignity.

- Reliable Quality
55gsm Polyethylene
breathable and anti-static material
- Sewn-on wrap-around tie for patient comfort and gown resilience
- Non transparent for patient confidence, privacy and dignity
- Overlocked edges for durability and gown strength
- Cost effective
- (colour may vary from previous stock)

CODE C-DIS1001



X-RAY EQUIPMENT



DEL RT100 PATIENT TABLE

The sturdy **RT100 four-way float top table** is suitable for use in a variety of general radiographic environments including small hospitals, clinics, ambulatory care centers and private offices. The table features a 700 lbs (318 kg) patient weight capacity and a stain resistant floating table top which glides freely on a precision roller bearing system. Float top movements are safely controlled by heavy-duty electric locks, activated by a conveniently placed foot treadle release. The 17" x 17" standard bucky tray can accommodate any cassette sized detector with ease. Its all steel structure and welded base provides maximum strength and reliability. In addition, the RT100's pedestal base can be used to store generator electronics making it a great space saver.

[RT100](#)

[Please contact us for pricing](#)

**CODE E-
GEN1041**



