

INFECTION CONTROL



DISPOSABLE PATIENT GOWN (BOX OF 100)

A simple change can have a large impact on patient dignity.

- Reliable Quality
55gsm Polyethylene
breathable and anti-static material
- Sewn-on wrap-around tie for patient comfort and gown resilience
- Non transparent for patient confidence, privacy and dignity
- Overlocked edges for durability and gown strength
- Cost effective
- (colour may vary from previous stock)

CODE C-DIS1001



X-RAY EQUIPMENT



DEL RT100 PATIENT TABLE

The sturdy
**RT100 four-way
float top table**

is suitable for use
in a variety of
general
radiographic
environments inc
luding small
hospitals, clinics,
ambulatory care
centers and
private offices.

The table
features a 700
lbs (318
kg) patient
weight capacity
and a stain
resistant floating
table top which
glides freely on a
precision roller
bearing system.

Float top
movements are
safely controlled
by heavy-duty
electric locks,
activated by a
conveniently
placed foot
treadle release.

The 17" x 17"
standard bucky
tray can
accommodate
any cassette
sized detector
with ease. Its all
steel structure
and welded base
provides
maximum
strength and
reliability. In
addition, the
RT100's pedestal
base can be used
to store
generator
electronics
making it a
great space
saver.

[RT100](#)

[Please contact us
for pricing](#)

**CODE E-
GEN1041**



