

POSITIONING AIDS



POSITIONING
SET OF MIXED
FOAMS (SET 2)

SET OF
MIXED
FOAMS

Variety of
shapes and
sizes
available, or
design your
own.

- 2x
FOAM0
2 -
Circle
- 2x
FOAM0
3 -
Small
Wedge
- 2x
FOAM0
4 - 15°
Wedge
- 1x
FOAM0
5 - 90°
Wedge
- 1x
FOAM0
6 -
Long 9
0°
Wedge
- 2x
FOAM0
7 -
Thin
Block

CODE C-
FOAM11-2



NECK INJECTION
POSITIONER CT

Nylon Neck
Injection
Positioning

Adjustable Strap.

For assisting
patients to keep
their shoulders
down during CT
contrast
injections. Strap
is positioned
under feet and
two loops held
tight to bring the
shoulders down.

CODE C-
PRO1038



SNAKE SANDBAG
1.3KG

100cm long x
6cm wide - 1.3kg

Double bagged
sandbag with a
hard wearing
Vinyl cover.

CODE C-SA1006



SNAKE SANDBAG
2.6KG

100cm long x
12cm wide -
2.6kg

Double bagged
sandbag with a
hard wearing
Vinyl cover.

CODE C-SA1007



**CUSTOM EASY
SLIDE, QUOTE
ON REQUEST**

Custom Patient
Easy Slide (non
friction

material for easy
sliding - made to
order)

CODE E-
ACC1035



**PATIENT EASY
SLIDE
LANDSCAPE**

Landscape Patient
Easy Slide (non
friction

material for easy
sliding).

CODE E-
ACC1036



**PATIENT EASY
SLIDE PORTRAIT**

Portrait Patient
Easy Slide (non
friction material for
easy sliding). Sold
Individually.

CODE E-
ACC1037

