

POSITIONING AIDS



POSITIONING  
SET OF MIXED  
FOAMS (SET 2)

SET OF  
MIXED  
FOAMS

Variety of  
shapes and  
sizes  
available, or  
design your  
own.

- 2x  
FOAM0  
2 -  
Circle
- 2x  
FOAM0  
3 -  
Small  
Wedge
- 2x  
FOAM0  
4 - 15°  
Wedge
- 1x  
FOAM0  
5 - 90°  
Wedge
- 1x  
FOAM0  
6 -  
Long 9  
0°  
Wedge
- 2x  
FOAM0  
7 -  
Thin  
Block

CODE C-  
FOAM11-2



NECK INJECTION  
POSITIONER CT

Nylon Neck  
Injection  
Positioning

Adjustable Strap.

For assisting  
patients to keep  
their shoulders  
down during CT  
contrast  
injections. Strap  
is positioned  
under feet and  
two loops held  
tight to bring the  
shoulders down.

CODE C-  
PRO1038



SNAKE SANDBAG  
1.3KG

100cm long x  
6cm wide - 1.3kg

Double bagged  
sandbag with a  
hard wearing  
Vinyl cover.

CODE C-SA1006



SNAKE SANDBAG  
2.6KG

100cm long x  
12cm wide -  
2.6kg

Double bagged  
sandbag with a  
hard wearing  
Vinyl cover.

CODE C-SA1007



**CUSTOM EASY  
SLIDE, QUOTE  
ON REQUEST**

Custom Patient  
Easy Slide (non  
friction

material for easy  
sliding - made to  
order)

**CODE E-  
ACC1035**



**PATIENT EASY  
SLIDE  
LANDSCAPE**

Landscape Patient  
Easy Slide (non  
friction

material for easy  
sliding).

**CODE E-  
ACC1036**



**PATIENT EASY  
SLIDE PORTRAIT**

Portrait Patient  
Easy Slide (non  
friction material for  
easy sliding). Sold  
Individually.

**CODE E-  
ACC1037**

