

POSITIONING AIDS



NECK INJECTION
POSITIONER CT

Nylon Neck
Injection
Positioning

Adjustable Strap.

For assisting
patients to keep
thier shoulders
down during CT
contrast
injections. Strap
is positioned
under feet and
two loops held
tight to bring the
shoulders down.

CODE C-
PRO1038



CUSTOM EASY
SLIDE, QUOTE
ON REQUEST

Custom Patient
Easy Slide (non
friction

material for easy
sliding - made to
order)

CODE E-
ACC1035



PATIENT EASY
SLIDE
LANDSCAPE

Landscape Patient
Easy Slide (non
friction

material for easy
sliding).

CODE E-
ACC1036



PATIENT EASY
SLIDE PORTRAIT

Portrait Patient
Easy Slide (non
friction material for
easy sliding). Sold
Individually.

CODE E-
ACC1037

