


POSITIONING AIDS



POSITIONING
SET OF MIXED
FOAMS (SET 1)

Foam wedges and blocks are a key component to creating high-quality radiographs.

- 2x FOAM0
2 - Circle:
100x50mm
- 2x FOAM0
3 - Small Wedge:
210x75x75mm
- 2x FOAM0
4 - 15° Wedge:
250x250x50mm
- 2x FOAM0
5 - 90° Wedge:
250x180x180mm
- 2x FOAM0
7 - Thin Block:
250x200x50mm
- 2x FOAM1
8 - Block:
200x100x100mm



POSITIONING
SET OF MIXED
FOAMS (SET 2)

SET OF
MIXED
FOAMS

Variety of shapes and sizes available, or design your own.

- 2x FOAM0
2 - Circle
- 2x FOAM0
3 - Small Wedge
- 2x FOAM0
4 - 15° Wedge
- 1x FOAM0
5 - 90° Wedge
- 1x FOAM0
6 - Long 90° Wedge
- 2x FOAM0
7 - Thin Block

CODE C-
FOAM11-2



0x100m
m
CODE C-
FOAM11-1

