

MARKERS



BEEKLEY R-SPOT® 113

"R" Spot
Metallic "R"
Orientation
marker
64 per box

CODE Beek113



**SPEE-D-MARK™
4MM CT SKIN
MARKER**

CT skin marker
Radiopaque
4mm pellet
50 per box


CODE SDM-CT2



**SPEE-D-
STRETCH 2.3MM
MAMMOGRAPHY
SKIN MARKER**

Nipple Marker
Super Stretchy
Radiopaque
2.3mm pellet


100 per box
CODE SDM-SS23



**SPEE-D-ANGLE™
ARROW
TOMOSYNTHESI
S
MAMMOGRAPHY
MARKER**

Palpable Mass
Marker
Radiopaque
Arrow Shape

100 per box
CODE SDM-TSD3



**SPEE-D-WING™
2MM
MAMMOGRAPHY
SKIN MARKER**

Nipple Marker
Radiopaque
2mm pellet
100 per box

CODE SDM-W20

