

MARKERS



**BEEKLEY  
TOMOSPOT®  
781**

**TOMOSPOT®  
FOR  
MARKING  
NIPPLES 781**

Distinctly identify and permanently document nipples, moles, scars, palpable masses, and areas of pain or concern on mammography images with BeekleySPOT® skin markers for mammography. The result is clear communication between patient, technologist and radiologist regarding areas of interest and more precise interpretation of the mammography.

**Ref: 781,  
(QTY:  
92/box)**

- 2.3mm non-metallic pellet
- pinch-free material with adhesive-free center



**BEEKLEY NEW  
TOMOSPOT®  
782**

**TOMOSPOT®  
FOR  
MARKING  
MOLES 782**

Provide certainty that densities on mammographic images highlighted by the O-shaped marker are moles and not areas of concern.

**Ref: 782,  
(QTY:  
54/box)**

- 1mm low-density see-through ring
- pinch-free material
- designed for use with 3D breast tomosynthesis equipment

**CODE C-  
BEEK782**



**BEEKLEY NEW  
TOMOSPOT®  
783**

**TOMOSPOT®  
FOR  
MARKING  
SCARS 783**

Distinctly identify and permanently document nipples, moles, scars, palpable masses, and areas of pain or concern on mammography images with BeekleySPOT® skin markers for mammography. The result is clear communication between patient, technologist and radiologist regarding areas of interest and more precise interpretation of the mammography.

**Ref: 783,  
(QTY:  
273cm cut  
measure  
roll/box)**

- 1mm diameter low-density see-through line
- pinch-free material



**BEEKLEY NEW  
TOMOSPOT®  
784**

**TOMOSPOT®  
FOR  
MARKING  
PALPABLE  
MASSES 784**

Distinctly identify and permanently document nipples, moles, scars, palpable masses, and areas of pain or concern on mammography images with BeekleySPOT® skin markers for mammography. The result is clear communication between patient, technologist and radiologist regarding areas of interest and more precise interpretation of the mammography.

**Ref: 784,  
(QTY:  
116/box)**

- low-density see-through triangle
- pinch-free material
- design



- design ed for use with 3D breast tomosynthesis equipment

CODE C-  
BEEK781

- design ed for use with 3D breast tomosynthesis equipment

CODE C-  
BEEK783

ed for use with 3D breast tomosynthesis equipment

CODE C-  
BEEK784



BEEKLEY NEW  
TOMOSPOT®  
785

**TomoSPOT®  
for marking  
non-palpable  
areas of  
concern / pain  
marker 785**

Distinctly identify and permanently document nipples, moles, scars, palpable masses, and areas of pain or concern on mammography images with BeeKleySPOT® skin markers for mammography. The result is clear communication between patient, technologist and radiologist regarding areas of interest and more precise interpretation of the mammography.

Ref: 785,  
(Qty: 58/box)

- low-density



see-through  
h  
square  
• pinch-free  
materi  
al  
• design  
ed for  
use  
with  
3D  
breast  
tomosy  
nthesis  
equip  
ment

**CODE C-**  
BEEK785

