

MARKERS



**BEEKLEY
TOMOSPOT®
781**

**TOMOSPOT®
FOR
MARKING
NIPPLES 781**

Distinctly identify and permanently document nipples, moles, scars, palpable masses, and areas of pain or concern on mammography images with BeekleySPOT® skin markers for mammography. The result is clear communication between patient, technologist and radiologist regarding areas of interest and more precise interpretation of the mammography.

**Ref: 781,
(QTY:
92/box)**

- 2.3mm non-metallic pellet
- pinch-free material with adhesive-free center



**BEEKLEY NEW
TOMOSPOT®
783**

**TOMOSPOT®
FOR
MARKING
SCARS 783**

Distinctly identify and permanently document nipples, moles, scars, palpable masses, and areas of pain or concern on mammography images with BeekleySPOT® skin markers for mammography. The result is clear communication between patient, technologist and radiologist regarding areas of interest and more precise interpretation of the mammography.

**Ref: 783,
(QTY:
273cm cut
measure
roll/box)**

- 1mm diameter low-density see-through line
- pinch-free material



- *design
ed for
use wit
h 3D
breast
tomosy
nthesis
equip
ment*

**CODE C-
BEEK781**

- *design
ed for
use
with
3D
breast
tomosy
nthesis
equip
ment*

**CODE C-
BEEK783**

