

DIAGNOSTIC EQUIPMENT



VET ENDOSCOPE

These VET Endoscope has special LED-AT-TIP technology in place of light guided fiber bundle:

1. **Super Bright LED-At-Tip:** Experience 30% more Brightness & 15% wider spread thanks to the new phosphor coating on the LED substrate, supporting high reliability over years of use.
2. **Fibreless LED-At-Tip:** Our patented **Fibreless LED-At-Tip technology** does away with the *fragile light guide glass fiber bundle* making the scopes very sturdy and environment friendly. Long life, superior light intensity, low maintenance



VIDEO VET ENDOSCOPE OVV-2H (L/I/S)

Video Vet Endoscope

Model: OVV-2H (L/I/S)

Super Bright LED-At-Tip: Experience 30% more Brightness & 15% wider spread thanks to the new Phosphor coating on the LED substrate, supporting high reliability over years of use.
HD Imaging Easy Manoeuverability

Remote Button for easy operation.

Ultra light weight, single step waterproof connector makes the endoscope more reliable & User-friendly.

CODE E-ESCOOVV-S2+US(1)-2

- nce & low replacement cost.
3. **Working life:** LED : 5000 Hrs and Xenon : 500 Hrs
 4. **Energy Consumption** is low in case of LED than compared to Xenon
 5. **Replacement Cost** is high due to shorter life span requiring more replacement in case of Xenon than compared to LED.
 6. **CMOS Sensor** provides unmatched image quality to explore smallest details.
 7. **Remote Button** for easy operation .
 8. **Ultra light weight,** single step waterproof connector makes the endoscope more reliable & User-friendly.

CODE E-ESCOOVV-S2+US(1)

