

## X-RAY EQUIPMENT



### RAYMASTER FULL SPINE INTEGRATED DIGITAL X-RAY CHIROPRACTIC SOLUTION

The **Raymaster Full Spine Integrated Digital X-Ray Chiropractic solution** is an economical and small footprint radiographic system designed specifically for chiropractic care. It combines flexibility, efficiency and affordability. The system will instantly improve patient workflow by eliminating the use of cassettes and developing conventional x-ray film. Instead, our Universal Chiropractic Solution are powered by the intuitive ATX Longtail full spine DR system giving technologists access to high resolution digital imaging and processing capabilities. The Longtail is the world's first monolithic full-spine flat-panel DR detector. No more long scan wait, no image stitching, no image gaps, and no plate overlap. The acquisition time is sub one second in most cases - so it will substantially improve efficiency.

#### Key benefits:



### RAYMASTER INTEGRATED DIGITAL X-RAY CHIROPRACTIC SOLUTION

The **Raymaster Integrated Digital X-Ray Chiropractic solution** is an economical and small footprint radiographic system designed specifically for chiropractic care. It combines flexibility, efficiency and affordability. The system will instantly improve patient workflow by eliminating the use of cassettes and developing conventional x-ray film. Instead, our Universal Chiropractic Solution are powered by the intuitive ATX DR system giving technologists access to high resolution digital imaging and processing capabilities. The acquisition time is sub one second in most cases - so it will substantially improve efficiency.

#### Key benefits of the Universal Upright Raymaster include:

- Mounts floor-to-wall or floor-to-ceiling
- Tube rotation enables precision



- Mounts floor-to-wall or floor-to-ceiling
- Tube rotation enables precision alignment of the tube for non-angled, angled and upright views.
- Convenient display gives indication of tube rotation
- CMP 200 X-ray generator series offers advanced high quality radiographic capabilities with unmatched performance and reliability resulting in the lowest cost of ownership over the product lifetime
- Longtail Digital Radiology. It is the first of its kind to produce long-bone and full-spine digital studies instantaneously with the highest geometric accuracy.

**Lead time:**  
Subject to stock availability, transit times and radiation licensing. Estimated shipping date to be provided once the purchase is confirmed

**CODE E-GEN1103**

- alignment of the tube for non-angled, angled and upright views.
- Convenient display gives indication of tube rotation
  - CMP 200 X-ray generator series offers advanced high quality radiographic capabilities with unmatched performance and reliability resulting in the lowest cost of ownership over the product lifetime.

**Lead time:**  
Subject to stock availability, transit times and radiation licensing. Estimated shipping date to be provided once the purchase is confirmed

**CODE EQU-GEN-10058**

