



# RADIATION PROTECTION



**LEAD GOWN  
PRICE  
CALCULATOR**

Choose your Size  
(XS) X Small: 18" / 46 cm Wide - Fits Chest 30-34" / 76.2-86.4cm  
(S) Small: 20" / 51 cm Wide - Fits Chest 34-38" / 86.4-96.5cm  
(M) Medium: 22" / 56 cm Wide - Fits Chest 38-42" / 96.5-106.7cm  
(L) Large: 24" / 61 cm Wide - Fits Chest 42-46" / 106.7-116.8cm

Regular lead =  
0.5mm

**CODE C-  
PRO1000**

## X-RAY EQUIPMENT



### DEL FLOOR MOUNTED TUBE STAND

The **FMT Floor mounted tube stand** features an advanced, 10.4" technologist-friendly touchscreen console that provides large text readouts and integrated operation of the tube stand. Ergonomically placed movement identification buttons provide optimal accuracy when setting up views. Take advantage of extensive longitudinal and vertical movements, as well as a 180° rotating tube column enabling positioning for cross-table lateral exams.

The FMT Floor Mounted Tube Stand comes in 3 different models:

[FMT18T](#)

[FMT18M](#)

[FMT](#)

**CODE E-GEN1032**



### DEL RT100 PATIENT TABLE

The sturdy **RT100 four-way float top table** is suitable for use in a variety of general radiographic environments including small hospitals, clinics, ambulatory care centers and private offices. The table features a 700 lbs (318 kg) patient weight capacity and a stain resistant floating table top which glides freely on a precision roller bearing system. Float top movements are safely controlled by heavy-duty electric locks, activated by a conveniently placed foot treadle release. The 17" x 17" standard bucky tray can accommodate any cassette sized detector with ease. Its all steel structure and welded base provides maximum strength and reliability. In addition, the RT100's pedestal base can be used to store generator electronics making it a great space saver.

[RT100](#)

[Please contact us for pricing](#)


**CODE E-GEN1041**



### RAYMASTER FULL SPINE INTEGRATED DIGITAL X-RAY CHIROPRACTIC SOLUTION

The **Raymaster Full Spine Integrated Digital X-Ray Chiropractic solution** is an economical and small footprint radiographic system designed specifically for chiropractic care. It combines flexibility, efficiency and affordability. The system will instantly improve patient workflow by eliminating the use of cassettes and developing conventional x-ray film. Instead, our Universal Chiropractic Solution are powered by the intuitive ATX Longtail full spine DR system giving technologists access to high resolution digital imaging and processing capabilities. The Longtail is the world's first monolithic full-spine flat-panel DR detector. No more long scan wait, no image stitching, no image gaps, and no plate overlap. The acquisition time is sub one second in most cases - so it will substantially improve efficiency.

**Key benefits:**



### ATX PORTABLE X-RAY GENERATOR 5KW

#### ATX Portable X-ray system

- Designed with innovation in mind, Ultra-light, high frequency x-ray system.
- Compact size and lightweight body
- Durable rugged exterior


Provided with Compliance testing as per local radiation health regulations.

#### Lead time:

Subject to stock availability, transit times and radiation licensing. Estimated shipping date to be provided once the purchase is confirmed

[Please contact us for pricing](#)

**CODE E-GEN1123**



- Mounts floor-to-wall or floor-to-ceiling
- Tube rotation enables precision alignment of the tube for non-angled, angled and upright views.
- Convenient display gives indication of tube rotation
- CMP 200 X-ray generator series offers advanced high quality radiographic capabilities with unmatched performance and reliability resulting in the lowest cost of ownership over the product lifetime
- Longtail Digital Radiology. It is the first of its kind to produce long-bone and full-spine digital studies instantaneously with the highest geometric accuracy.

**Lead time:**  
Subject to stock availability, transit times and radiation licensing. Estimated shipping date to be provided once the purchase is confirmed

**CODE E-**  
GEN1103



**ATX PORTABLE  
X-RAY  
GENERATOR  
2.4KW**

**ATX Portable X-ray system**

- Designed with innovation in mind, Ultra-light, high frequency x-ray system.
- Compact size and lightweight body
- Durable rugged exterior

Provided with Compliance testing as per local radiation health regulations.

**Lead time:**  
Subject to stock availability, transit times and radiation licensing. Estimated shipping date to be provided once the purchase is confirmed

[Please contact us for pricing](#)

**CODE E-  
GEN1124**



**ATX PORTABLE  
X-RAY  
GENERATOR  
1.8KW**

**ATX Portable X-ray system**

- Designed with innovation in mind, Ultra-light, high frequency x-ray system.
- Compact size and lightweight body
- Durable rugged exterior

Provided with Compliance testing as per local radiation health regulations.

**Lead time:**  
Subject to stock availability, transit times and radiation licensing. Estimated shipping date to be provided once the purchase is confirmed

[Please contact us for pricing](#)

**CODE E-  
GEN1125**



**ATX SMART  
PORTABLE X-  
RAY GENERATOR**

**ATX Smart-portable X-ray system**

- Embedded tablet touch screen console
- Designed with innovation in mind, Ultra-light, high frequency x-ray system.
- Compact size and lightweight body
- Durable rugged exterior
- Auto calibration function
- Body part APR function (over 960 APR)
- DAP display function
- Also available with a traditional button keypad control

Compliance testing as per local radiation health regulations provided

**Lead time:**  
Subject to stock availability, transit times and radiation licensing. Estimated shipping date to be provided once the purchase is confirmed

**CODE EQU-  
GEN-10020**



**RAYMASTER  
INTEGRATED  
DIGITAL X-RAY  
CHIROPRACTIC  
SOLUTION**

The **Raymaster Integrated Digital X-Ray Chiropractic solution** is an economical and small footprint radiographic system designed specifically for chiropractic care. It combines flexibility, efficiency and affordability. The system will instantly improve patient workflow by eliminating the use of cassettes and developing conventional x-ray film. Instead, our Universal Chiropractic Solution are powered by the intuitive ATX DR system giving technologists access to high resolution digital imaging and processing capabilities. The acquisition time is sub one second in most cases - so it will substantially improve efficiency.

**Key benefits of the Universal Upright Raymaster include:**

- Mounts floor-to-wall or floor-to-ceiling
- Tube rotation enables precision alignment of the tube for non-angled, angled and upright views.
- Convenient display gives

•

**Lead time:**  
Subject to stock  
availability,  
transit times and  
radiation  
licensing.  
Estimated  
shipping date to  
be provided once  
the purchase is  
confirmed

**Shop online at [www.atxmedical.com.au](http://www.atxmedical.com.au)**

## X-RAY STATIONERY



**X-RAY  
ENVELOPES  
444X368MM  
WHITE FOR  
35X43 FILM**

Box of 250

**CODE C-**  
GEN1004