

## MARKERS



### BEEKLEY X-SPOT® 101

#### BEEKLEY X-SPOT 101 BLUE

Ref: 101 Blue,  
(QTY: 150/box)

*X-SPOT 101: Identify and differentiate specific points of reference in conventional simulation treatment planning, chest, and trauma imaging.*

- 1.5mm lead-free pellet
- easy Pick-Up tab

**CODE C-**  
BEEK101



### BEEKLEY Y-SPOT® 102

#### BEEKLEY Y-SPOT 102 Green

Ref: 102 Green  
(QTY: 132/box)

*Identify and differentiate specific points of reference in conventional simulation treatment planning, chest, and trauma imaging.*

- 2.0mm lead-free pellet
- easy Pick-Up tab

**CODE C-**  
BEEK102



### BEEKLEY D-SPOT® ARROW 104

#### BEEKLEY D-SPOT® ARROW 104

Ref: 104 (QTY: 99/box)

- lead-free arrow with easy Pick-Up tab
- for marking nodes and outer canthus

**CODE C-**  
BEEK104



### BEEKLEY T-SPOT® LINE 111

#### BEEKLEY T-SPOT® LINE 111

Ref: 111, (QTY: 275cm cut to measure roll/box)

- .8mm diameter metallic line.
- for marking scars, nodes, field borders, and match lines

**CODE C-**  
BEEK111



**BEEKLEY  
GUIDELINES® CT  
BIOPSY GRID  
117**

**Guidelines® CT  
Biopsy  
Grid 117**  
**Ref: 117,  
(QTY:  
20/box)**

- 4" x 5" biopsy grid
- 7 non-metallic radiopaque lines
- open pane design for easy mark placement
- ideal for CT-guided biopsies of the neck and other small areas

**CODE C-  
BEEK117**



**BEEKLEY CT-  
SPOT® 119**

**BEEKLEY CT-  
SPOT® 119**

**Ref: 119, (QTY:  
85/box)**

*Ideal for marking 3-point set-ups and isocenters, these unique markers image brightly during simulation without streaking or artifact.*

- 2.3mm non-metallic pellet
- for 3-points set-ups and isocenters

**CODE C-  
BEEK119**




**BEEKLEY CT-  
SPOT® 120**

*Ideal for marking 3-point set-ups and isocenters, these unique markers image brightly during simulation without streaking or artifact.*

**BEEKLEY CT-  
SPOT® 120,  
Ref: 120 (QTY:  
85/box)**

- 4.0mm non-metallic pellet
- for 3-point set-ups and isocenters

**CODE C-  
BEEK120**



**BEEKLEY  
PINPOINT®128**

**PinPoint® for  
Image  
Registration  
128**

**Ref: 128, (QTY:  
30 / box)**

- ideal for image registration performed under CT and MRI
- conical design with a 1.27mm diameter center hole

**CODE C-  
BEEK128**



**MR-SPOT  
PACKETS™ 185**

**MR-SPOT  
Packets™ 185**

**Ref: 185, (QTY:  
40/box)**

- 1.75cm  
Radiance  
® filled  
packet
- ideal for  
extremities,  
soft  
tissue  
masses,  
and  
breast  
MRI
- designed  
for  
optimum  
patient  
comfort

**CODE C-  
BEEK185**



**BEEKLEY  
LOCALIZERS™  
210**

**Localizer  
Localizers™  
on  
Markers  
210**

**Ref: 210,  
(QTY:  
40/box)**

- 2.3mm  
low-  
density  
pellet
- ideal  
for  
identifying  
soft  
tissue  
masses  
on the  
neck

**CODE C-  
BEEK210**



**BEEKLEY  
LOCALIZERS™  
211**

**LOCALIZERS  
™  
LOCALIZATI  
ON  
MARKERS  
211**

**Ref: 211,  
(QTY:  
40/box)**

- 4.0mm  
low-  
density  
pellet
- ideal  
for  
identifying  
areas  
of  
interest  
on  
abdominal  
cases  
or with  
larger  
patients

**CODE C-  
BEEK211**



**BEEKLEY  
POINTGUARDS®  
215**

**PointGuards®  
215**

**Ref: 215, (QTY:  
75/box)**

- 1.5cm x  
5.0cm  
linear  
adhesive  
covers
- for  
longer,  
more  
linear  
treatment  
planning  
marks

**CODE C-  
BEEK215**



**BEEKLEY  
GUIDELINES® CT  
BIOPSY GRID  
217**

**Guidelines® CT  
Biopsy  
Grid 217**

**Ref: 217,  
(QTY: 15 /  
box)**

- 7" x 7.5" biopsy grid
- 12 non-metallic radiopaque lines
- open pane design for easy mark placement
- ideal for CT-guided biopsies of abdomen and other large areas

**CODE C-  
BEEK217**



**BEEKLEY  
ACCUGRID®300  
4.65"**

**AccuGrid®  
Specimen  
Localizing Grid  
System 300**

**Ref: 300, (QTY:  
12 / box)**

- 4.65" radiolucent grid inside a sealable surgical bag
- for the safe imaging, transport, and processing of surgical breast specimens

**CODE C-  
BEEK300**



**BEEKLEY  
ACCUGRID®302  
5.65"**

**AccuGrid®  
Specimen  
Localizing Grid  
System 302**

**Ref: 302, (QTY:  
10 / box)**

- 5.65" radiolucent grid inside a sealable surgical bag
- for the safe imaging, transport, and processing of surgical breast specimens
- expanded grid size for use with larger specimens

**CODE C-  
BEEK302**



**BEEKLEY  
CRADLES®309**

**Cradles®  
Needle  
Localization  
Wire  
Protectors 309**

**Ref: 309, (QTY:  
12 / box)**

- 12" x 3" translucent cover with foam adhesive border
- disposable, clear covers for protecting and stabilizing the needle localization wire

**CODE C-  
BEEK309**



BEEKLEY BELLA  
BLANKETS®310

**Bella  
Blankets  
®  
protective  
coverlets  
310**

Provide  
sanitary  
protection  
for patients  
and help  
immobilize  
hard to  
position  
breasts with  
Bella  
Blankets. Im-  
prove  
Clinical  
Outcomes,  
Increase  
Patient  
Satisfaction,  
Save Money  
and Increase  
Revenues.

**Ref: 310,  
(QTY:  
50/box)**

- 22.9cm x 27.9cm
- intended for use on 18cm x 24cm recept or plate

**CODE C-  
BEEK310**




BEEKLEY  
POINTGUARDS®  
315

**PointGuards® 315**

**Ref: 315,  
(QTY: 100 /  
box)**

- 2.5cm circular adhesive covers
- for smaller treatment planning marks
- see-through cover for quicker and more precise placement

**CODE C-  
BEEK315**



BEEKLEY  
POINTGUARDS®  
318

**PointGuards®  
318**

**Ref: 318, (QTY:  
85/box)**

- 2.5cm circular adhesive covers
- crosshair design
- for hard to find marks

**CODE C-  
BEEK318**



BEEKLEY  
COLDSPOT® 500

**ColdSPOT®  
500**

**Sold  
Individually**

- 2.75" x 3.0" gel pack
- self-adhesive side tabs
- for targeted and effective cold treatment following minimally invasive breast biopsy procedures

**CODE C-  
BEEK500**



**BEEKLEY X-SPOT® 601 FLORALS®**

**X-SPOT® 601  
Florals®  
101  
Blue**

**Ref: 601  
Florals®  
Nipple  
Marker,  
(QTY:  
150/box)  
Available  
in blue  
(Ref:101)**

- 1.5mm lead-free pellet
- easy Pik-UP tab
- compatible with both digital and analog equipment

**CODE C-BEEK601**



**BEEKLEY S-SPOT® 607 FLORALS®**

**S-SPOT® CT  
Simulation  
Marker 607  
Florals®**

Help define architectural distortion at the biopsy site by delineating the exact location of incision.

**Ref: 607, (QTY:  
363cm cut to  
measure  
roll/box)**

- 33mm diameter radiopaque line
- burnout-resistant brighter image
- for marking scars and field borders

**CODE C-BEEK607**



**BEEKLEY A-SPOT® 651 FLORALS®**

**A-SPOT® Light  
Image®  
palpable mass  
marker 651  
Florals®**

**Ref: 651 Floral  
s®, (QTY:120  
/ box)**

- uniform low-density see-through triangle
- highlights without obscuring
- compatible with both digital and analog equipment

**CODE C-BEEK651**



**BEEKLEY O-SPOT® LIGHT IMAGE® 652 FLORALS**

**O-SPOT®  
Light  
Image®  
652  
Florals**

**Ref: 652  
(QTY:58/  
box)**

- 1mm uniform low-density see-through ring
- highlights without obscuring
- compatible with both digital and analog equipment

**CODE C-BEEK652**



**BEEKLEY S-SPOT® 677 FLORALS®**

**S-SPOT® Light Image® Soft 'n Stretchy® 677 Florals®**

**Ref: 677 Florals®, (QTY: 280cm cut to measure roll / box)**

- 1mm diameter uniform low-density see-through line
- pinch-free material
- compatible with both digital and analog equipment

**CODE C-BEEK677**



**BEEKLEY V-SPOT®703 FLORALS®**

**V-SPOT® Designed for Digital®**

**Ref: 703 Florals®, (QTY: 130/box)**

- 2.3mm non-metallic pellet
- easy Pik-UP tab
- designed for use with the most sensitive digital equipment

**CODE C-BEEK703**



**BEEKLEY TOMOSPOT® 781**

**TOMOSPOT® FOR MARKING NIPPLES 781**

Distinctly identify and permanently document nipples, moles, scars, palpable masses, and areas of pain or concern on mammography images with BEEKLEYSPOT® skin markers for mammography. The result is clear communication between patient, technologist and radiologist regarding areas of interest and more precise interpretation of the mammography.

**Ref: 781, (QTY: 92/box)**

- 2.3mm non-metallic pellet
- pinch-free material with adhesive center
- designed for use with 3D breast tomosynthesis equipment

**CODE C-BEEK781**



**BEEKLEY NEW TOMOSPOT® 782**

**TOMOSPOT® FOR MARKING MOLES 782**

Provide certainty that densities on mammographic images highlighted by the O-shaped marker are moles and not areas of concern.

**Ref: 782, (QTY: 54/box)**

- 1mm low-density see-through ring
- pinch-free material
- designed for use with 3D breast tomosynthesis equipment

**CODE C-BEEK782**





**BEEKLEY NEW  
TOMOSPOT®  
783**

**TOMOSPOT®  
FOR  
MARKING  
SCARS 783**

Distinctly identify and permanently document nipples, moles, scars, palpable masses, and areas of pain or concern on mammography images with BeekleySPOT® skin markers for mammography. The result is clear communication between patient, technologist and radiologist regarding areas of interest and more precise interpretation of the mammography.

**Ref: 783,  
(QTY:  
273cm cut  
measure  
roll/box)**

- 1mm diameter low-density see-through line
- pinch-free material
- designed for use with 3D



**BEEKLEY NEW  
TOMOSPOT®  
784**

**TOMOSPOT®  
FOR  
MARKING  
PALPABLE  
MASSES 784**

Distinctly identify and permanently document nipples, moles, scars, palpable masses, and areas of pain or concern on mammography images with BeekleySPOT® skin markers for mammography. The result is clear communication between patient, technologist and radiologist regarding areas of interest and more precise interpretation of the mammography.

**Ref: 784,  
(QTY:  
116/box)**

- low-density see-through triangle
- pinch-free material
- designed for use with 3D



**BEEKLEY NEW  
TOMOSPOT®  
785**

**TomoSPOT®  
for marking  
non-palpable  
areas of  
concern / pain  
marker 785**

Distinctly identify and permanently document nipples, moles, scars, palpable masses, and areas of pain or concern on mammography images with BeekleySPOT® skin markers for mammography. The result is clear communication between patient, technologist and radiologist regarding areas of interest and more precise interpretation of the mammography.

**Ref: 785,  
(QTY:  
58/box)**

- low-density see-through square
- pinch-free material
- designed for use with 3D breast tomosy



**BEEKLEY N-  
SPOT® 790  
FLORALS®**

**N-SPOT®  
Designed  
for  
Digital®  
Soft n  
Stretchy  
®**

**Ref: 790  
Florals®,  
(QTY:  
94/box)**

- 2.3mm non-metallic pellet
- pinch-free material with adhesive-free center
- designed for use with the most sensitive digital equipment

**CODE C-  
BEEK790**





breast  
tomosy  
nthesis  
equip  
ment

CODE C-  
BEEK783

breast  
tomosy  
nthesis  
equip  
ment

CODE C-  
BEEK784

nthesis  
equip  
ment

CODE C-  
BEEK785



BEEKLEY O-  
SPOT®791  
FLORALS®

O-SPOT®  
Designed  
for  
Digital®  
Soft n  
Stretchy  
®791  
Florals®

Provide  
certainty  
that  
densities on  
mammograp  
hic images  
highlighted  
by the O-  
shaped  
marker are  
moles and  
not areas of  
concern.

Ref: 791  
(QTY:  
55/box)

- 1mm  
unifor  
m  
ultra-  
sheer  
see-  
throug  
h ring
- pinch-  
free  
materi  
al
- design  
ed for  
use  
with  
the  
most  
sensiti  
ve  
digital  
equip  
ment

CODE C-  
BEEK791



