

MARKERS




**AKUCAL IMAGE
SCALING DEVICE
(25MM)
(COMPLETE)**

25mm image
scaling device

COMPLETE with
all parts
including-

- Suction cup base for Horizontal or Vertical use
- Flexible arm to position marker at anatomy level
- Aluminium tip avoids breakage from wear and tear
- Numerical size indicators visible in every image (25mm, 25.4mm, 30mm)
- Surgical Grade stainless steel marker
- Lever to engage or disengage suction

**CODE E-
ACC1029-25**



**AKUCAL
SCALING TOOL (25MM
REPLACEMENT
MARKER)**

25mm

Replacement
stainless steel
marker for Image
Scaling Tool

EACH

**CODE E-
ACC1030-25**



**AKUCAL IMAGE
SCALING DEVICE
(30MM)
(COMPLETE)**

25mm image
scaling device

COMPLETE with
all parts
including-

- Suction cup base for Horizontal or Vertical use
- Flexible arm to position marker at anatomy level
- Aluminium tip avoids breakage from wear and tear
- Numerical size indicators visible in every image (25mm, 25.4mm, 30mm)
- Surgical Grade stainless steel marker
- Lever to engage or disengage suction

**CODE E-
ACC1031-30**



**AKUCAL
SCALING TOOL (30MM
REPLACEMENT
MARKER)**

30mm

Replacement
stainless steel
marker for Image
Scaling Tool

EACH

**CODE E-
ACC1032-30**

